NADAL[®] Malaria Pf/Pan Ag 4 Species Test (test cassette)

REF 472038N-25

de Gebrauchsanweisung	2	pt Instruções de Utilização	26
en Instructions for use	6	sv Användarinstruktioner	30
fr Instructions d'utilisation	10	nl Gebruiksaanwijzing	34
es Instrucciones de uso	14	Symbols	39
it Istruzioni per l'uso	18	Our Teams	40
pl Sposób użycia	22		

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Rev. 0, 2023-11-13 OM/SDe

1. Intended Use

The NADAL® Malaria Pf/Pan Ag 4 Species Test is a lateral flow chromatographic immunoassay for the simultaneous, qualitative detection and differentiation of HRP II (histidine-rich protein II) specific to *Plasmodium falciparum* (*P. falciparum*) and pLDH (plasmodium lactate dehydrogenase) specific to *Plasmodium* species (*P. falciparum*, *P. vivax*, *P. malariae*, *P. ovale*) in human whole blood specimens. The test is intended for use as an aid in the diagnosis of malaria infection and is designed for professional use only. Any specimen found to be reactive with the NADAL® Malaria Pf/Pan Ag 4 Species Test should be confirmed with a supplementary assay, e. g. the microscopic examination of a thin blood smear.

2. Introduction and Clinical Significance

Malaria is a mosquito-borne, haemolytic, febrile illness. It infects over 200 million people and kills more than 400,000 people annually.¹ Malaria is caused by four species of *Plasmodium*: *P. falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. All these plasmodia infect and destroy human erythrocytes, causing such symptoms as chills, fever, anaemia and splenomegaly. *P. falciparum* causes more severe disease than the other *Plasmodium* species and accounts for most malaria-related deaths. *P. falciparum* and *P. vivax* are the most common pathogens; there is, however, a considerable geographic variation in the distribution of *Plasmodium* species.²

Traditionally, malaria is diagnosed by the detection of organisms on Giemsa-stained thick smears of peripheral blood, and different *Plasmodium* species are differentiated by their presence in infected erythrocytes.² The technique provides an accurate and reliable diagnosis, but only when performed by skilled microscopists using defined protocols,³ which presents major hindrances for remote and poor areas of the world.

The NADAL® Malaria Pf/Pan Ag 4 Species Test enables the simultaneous, qualitative detection and differentiation of HRP II specific to *P. falciparum* and/or pLDH produced by all *Plasmodium* species. The test utilises colloidal gold conjugates to selectively detect HRP II and pLDH in human whole blood samples.

In 2010, researchers working in the Peruvian Amazon region identified patients infected with *P. falciparum* strains that carried HRP II and HRP III gene deletions, rendering them undetectable by HRP II-based rapid tests. The frequency and global distribution of this phenomenon is not yet fully understood, but in a limited number of countries, the relative incidence of these deletion mutants has been found to be high enough to question the usefulness of rapid tests detecting only HRP II.

3. Test Principle

The NADAL® Malaria Pf/Pan Ag 4 Species Test is a lateral flow chromatographic immunoassay for the simultaneous, qualitative detection and differentiation of HRP II specific to *P. falciparum* and pLDH specific to *Plasmodium* species (*P. falciparum*, *P. vivax*, *P. malariae*, *P. ovale*) in human whole blood specimens. The membrane is pre-coated with anti-HRP II antibodies and anti-pLDH antibodies.

During testing, the whole blood sample reacts with the anti-HRP II antibodies and anti-pLDH antibodies-dye conjugates, which are pre-coated on the internal test strip. The mixture then migrates along the membrane by capillary action and reacts with anti-HRP II antibodies in the test line region 'Pf' and with anti-pLDH antibodies in the test line region 'Pan'. If the sample contains HRP II or/and pLDH, a coloured line will form in the test line region 'Pf' or 'Pan', or a coloured line will form in both the 'Pf' test line region and the 'Pan' test line region.

The absence of a coloured line in the test line region 'Pf' and/or 'Pan' indicates that the sample does not contain HRP II and/or pLDH.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of sample has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 25 NADAL® Malaria Pf/Pan Ag 4 Species test cassettes
- 25 disposable pipettes (10 µL)
- 1 buffer (7 mL)
- 1 package insert

5. Additional Materials Required

- Timer
- Sample collection container
- Lancets (for fingerstick whole blood only)
- Alcohol pads
- Micropipette (if necessary)

6. Storage & Stability

NADAL® Malaria Pf/Pan Ag 4 Species test kits should be stored at 2-30°C and used by the expiry date indicated on the packaging. Test cassettes should remain in sealed pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date. Care should be taken to protect components of the test kit from contamination. Do not use the test if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiry date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of samples by using a new sample collection container for each sample obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where samples and test kits are handled.

- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when samples are being assayed.
- Handle all samples as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of samples.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

The NADAL® Malaria Pf/Pan Ag 4 Species Test can be performed using whole blood (from venipuncture or fingerstick).

To collect fingerstick whole blood samples:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood samples

Containers containing anticoagulants, such as EDTA, sodium citrate, potassium oxalate or sodium heparin should be used for the preparation of venous whole blood samples.

Testing should be performed immediately after sample collection. Do not leave samples at room temperature for prolonged periods of time.

If the test is to be run within 3 days of sample collection, whole blood collected by venipuncture should be stored at 2-8°C. For long-term storage, samples should be kept at -20°C.

Bring samples to room temperature prior to testing. Frozen samples should be completely thawed and mixed well prior to testing. Samples should not be frozen and thawed repeatedly.

If samples are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

9. Test Procedure


Bring tests, samples, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch.

Label the test cassette with the patient or control identification.

2. Place the test cassette on a clean and level surface.


3. Holding the pipette vertically, draw the sample up to the first widening (approximately 10 µL) and add it to the upper well (W1) of the test cassette, then add 3 full drops of buffer to the lower well (W2). Avoid air bubbles.

4. Start the timer.

5. After exactly 5 minutes, add 1 full drop of buffer to the upper well (W1).

6. Wait for the coloured line(s) to appear. Read the test result after a further 10 minutes. Do not interpret the result after more than 20 minutes in total.

10. Result Interpretation


Positive for *P. falciparum* infection

Two coloured lines appear in the result area. One coloured line appears in the control line region (C) and the other coloured line appears in the test line region 'Pf'. The test result indicates the presence of HRP II and therefore, a *P. falciparum* infection.

Positive for non-*P. falciparum* infection (*P. vivax*, *P. malariae* or *P. ovale*)*

Two coloured lines appear in the result area. One coloured line appears in the control line region (C) and the other coloured line appears in the test line region 'Pan'. The test result indicates the presence of pLDH and therefore, a *P. vivax*, *P. malariae* or for *P. ovale* infection.

* in infections with *P. falciparum* HRP II deletion mutants the line in the test line region 'Pf' does not appear. Therefore, in regions endemic for HRP II deletion mutants, this kind of result could also indicate an infection with the *P. falciparum* deletion mutant alone or a mixed infection. A differentiation between infection types is not possible.

Positive for *P. falciparum*/mixed *Plasmodium* infection

Three coloured lines appear in the result area. One line appears in the control line region (C) and one coloured

line appears in each test line region 'Pf' and 'Pan'. The test result indicates the presence of both HRP II and pLDH, and therefore, either a single infection with *P. falciparum* or infection with several *Plasmodium*-species (*P. vivax*, *P. malariae* and/or *P. ovale*), including *P. falciparum*.

Note: Samples with positive results should be confirmed with alternative testing method(s) and clinical findings before a positive determination is made.

Negative

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region 'Pf' and 'Pan', which indicates that no *Plasmodium* antigens have been detected.

Invalid

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note:

The colour intensity in the test line region may vary depending on the concentration of the analyte present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of the analyte in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Malaria Pf/Pan Ag 4 Species Test is for professional *in-vitro* diagnostic use only. The test should only be used for the qualitative detection of *P. falciparum*, *P. vivax*, *P. malariae* or *P. ovale* antigens in human whole blood samples. Neither the quantitative value nor the rate of increase in the concentration of *Plasmodium* parasites can be determined with this qualitative test.

- As described in section 10 'Result Interpretation', infection with *P. falciparum* HRP II deletion mutants cannot be distinguished from infection with non-*P. falciparum* species. Therefore, in areas where such mutants are endemic, additional testing must be done, especially if a clear differentiation is required for patient management decisions.
- The NADAL® Malaria Pf/Pan Ag 4 Species Test only indicates the presence of *P. falciparum*, *P. vivax*, *P. malariae* or *P. ovale* antigens in the sample and should not be used as the sole criterion for the diagnosis of a malaria infection.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a malaria infection.

13. Performance Characteristics

Analytical sensitivity

The NADAL® Malaria Pf/Pan Ag 4 Species Test for the detection of *P. falciparum* infection and non-*P. falciparum* infection according to the parasite level:

	No. of parasites/ μ L of blood	Microscopy	NADAL® Malaria Pf/Pan Ag 4 Species Test	Sensitivity
<i>P. falciparum</i>	1-100	12	11	91.7%
	101-200	15	15	100%
	201-500	10	10	100%
	500-1000	9	9	100%
	>1000	7	7	100%
	Total	53	52	98.1%
Non- <i>P. falciparum</i>	1-100	38	37	97.4%
	101-200	43	43	100%
	201-500	37	37	100%
	500-1000	28	28	100%
	>1000	12	12	100%
	Total	158	157	99.4%

Clinical sensitivity and specificity

Method		Reference method (microscopy)			
		Positive		Negative	Total
NADAL® Malaria Pf/Pan Ag 4 Species Test	Positive	Pv	Pf		
				157	0
		0	52	3	55
	Negative	1	1	324	326
	Total	158	53	327	538

Note: The comparison study for the 'Pan' line was carried out using blood specimens positive only for *P. vivax*.

Relative sensitivity: $(157+52)/(158+53)=99.1\%$ (96.6%-99.7%)

Relative specificity: $324/327=99.1\%$ (97.3%-99.7%)

Overall agreement: $(157+52+324)/(158+53+327)=99.1\%$
(97.8%-99.6%)

*95% Confidence interval

Precision

Reproducibility was established by testing triplicates of 6 specimens. Repeatability was established within the reproducibility study. Testing was performed by 3 operators, using 3 independent NADAL® Malaria Pf/Pan Ag 4 Species test lots at 3 testing sites on 5 separate days.

The NADAL® Malaria Pf/Pan Ag 4 Species Test demonstrated acceptable repeatability and reproducibility when performed by trained professionals. The negative and positive values were correctly identified >99% of the time.

Interfering substances

Specimens were spiked with the following potentially interfering substances and evaluated using the NADAL® Malaria Pf/Pan Ag 4 Species Test.

Potentially interfering substances	Concentration
Haemoglobin	2 g/L
Bilirubin (conjugated)	342 µmol/L
Bilirubin (unconjugated)	342 µmol/L
Triglycerides	37 mmol/L
Protein (total)	120 g/L
Albumin	60 g/L
γ globulin	60 g/L
Human anti-mouse antibody (HAMA)	780.73 IU/mL
Rheumatoid factor	1035 IU/mL
Metronidazole	701 µmol/L
Quinine	148 µmol/L
Rifampicin	78.1 µmol/L
Aspirin	4.34 mmol/L
Paracetamol	199 µmol/L
Ibuprofen	2425 µmol/L
Caffeine	308 µmol/L
Ethanol	86.8 mmol/L

None of the substances interfered with the assay at the concentrations tested.

Cross-reactivity

Clinical specimens were evaluated for cross-reactivity with unrelated infections using the NADAL® Malaria Pf/Pan Ag 4 Species Test.

Unrelated infections	No. of specimens
anti-HIV+	10
HBsAg+	20
anti-HCV+	10
anti-HAV+	20
anti-Rubella+	20
anti-CMV+	20
anti-HSV-1+	20
anti-HSV-2+	20
Typhoid fever+	20
Yellow fever+	5
Dengue fever+	5

Unrelated infections	No. of specimens
anti-Influenza A+	20
anti-Influenza B+	10
anti-Toxoplasma gondii+	10
Syphilis+	20
Tuberculosis+	10

No cross-reactivity was observed with the aforementioned unrelated infections when tested using the NADAL® Malaria Pf/Pan Ag 4 Species Test.

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