NADAL[®] Dengue NS1 Ag + IgG/IgM Test (test cassette)

REF 532919N-10

de	Gebrauchsanweisung	2	cs	Návod k použití	37
en	Instructions for use	7	fi	Käyttöohje	42
fr	Instructions d'utilisation	12		Symbols	47
es	Instrucciones de uso	17		Our Teams	48
it	Istruzioni per l'uso	22			
pl	Sposób użycia	27			
pt	Instruções de Utilização	32			

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15. Referenzen

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German).* 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

Rev. 1, 2023-01-19 OM

1. Intended Use

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of dengue virus NS1 antigens as well as anti-dengue virus IgG and IgM in human whole blood, serum or plasma specimens (see section 12 'Limitations'). The test is intended for use as an aid in the diagnosis of primary and secondary dengue virus infections in patients with febrile disease of unknown origin mainly associated with tropical regions, who are showing symptoms such as sudden fever with intense headaches and limb pain (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Dengue NS1 Ag + IgG/IgM Test is designed for professional use only.

2. Introduction and Clinical Significance

Dengue virus is a flavivirus transmitted primarily by tiger mosquitoes *Aedes aegypti* and *Aedes albopictus*. It is widely distributed throughout the tropical and subtropical areas of the world and causes up to 100 million infections annually. Classic dengue infection is characterised by a sudden onset of fever, intense headache, myalgia, arthralgia and a rash.

Primary dengue infection causes IgM to increase to a detectable level approximately 5 days after the onset of fever. IgM generally persists for 30 to 90 days. Most dengue patients in endemic regions have secondary infections, resulting in high levels of specific IgG prior to, or simultaneous with, IgM response. Therefore, the detection of specific anti-dengue IgM and IgG can also help to distinguish between primary and secondary infections.

NS1 is one of 7 dengue virus non-structural proteins which are thought to be involved in viral replication. NS1 exists as a monomer in its immature form but is rapidly processed in the endoplasmic reticulum to form a stable dimer. A small amount of NS1 remains associated with intracellular organelles, where it is thought to be involved in viral replication. The rest of NS1 is found to be either associated with the plasma membrane or secreted as a soluble hexadimer. NS1 is essential for viral viability but its precise biological function is unknown. Antibodies produced in response to NS1 in viral infections can cross-react with cell surface antigens on epithelial cells and platelets. This may have implications in the development of dengue haemorrhagic fever.

3. Test Principle

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of dengue virus NS1 antigens, as well as anti-dengue virus IgG and IgM in human whole blood, serum or plasma specimens.

The NADAL® Dengue NS1 Ag + IgG/IgM Test contains two internal test strips (left side of the test cassette: Dengue NS1 Ag Test; right side of the test cassette: Dengue IgG/IgM Test).

Anti-dengue virus NS1 antibodies are immobilised in the test line region 'T' of the membrane. Anti-dengue virus IgM and anti-dengue virus IgG are immobilised in the test line regions of the membrane labelled 'IgM' and 'IgG' respectively.

During the test, dengue virus NS1 antigens, if present in the specimen, bind to anti-dengue virus NS1 antibodies, which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. If specific anti-dengue virus IgM or anti-dengue virus IgG are present in the specimen, they bind to the dengue virus antigens conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-dengue virus NS1 antibodies as well as anti-dengue virus IgM and anti-dengue virus IgG in the test line regions 'T' as well as 'IgM' and 'IgG' respectively. Excess coloured particles are captured in each control line region 'C'. The presence of a coloured line in the test line region 'T' and/or 'IgM' and/or 'IgG' indicates a positive result. The absence of a coloured line in the test line region 'T' and/or 'IgM' and/or 'IgG' indicates a negative result.

The formation of a coloured line in each control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® Dengue NS1 Ag + IgG/IgM test cassettes, incl. disposable pipettes (40 µL) for dengue NS1 Ag
- 10 capillary tubes (10 µL) for dengue IgG/IgM
- 1 buffer (3 mL)*
- 1 package insert

*containing the following preservative: sodium azide: <0.1%.

No hazard labelling for sodium azide is required according to Regulation (EC) No 1272/2008 CLP. Concentrations are below the exemption threshold of 0.1%.

5. Additional Materials Required

- Specimen collection containers (appropriate for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.

- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

The NADAL® Dengue NS1 Ag + IgG/IgM Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood specimens

Containers containing anticoagulants, such as EDTA, sodium citrate or heparin sodium should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 24 hours of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 5 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure


Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
2. Place the test cassette on a clean and level surface.

3. a) Dengue virus NS1 Ag


- Holding a pipette vertically, transfer 1 drop (approximately 30-45 µL) of the whole blood, serum or plasma specimen to the specimen well (S1) of the test cassette.
- Holding the buffer bottle vertically, add 1 drop (approximately 30-40 µL) of buffer to the specimen well (S1).

Avoid trapping air bubbles in the specimen well (S1) and do not add any solution to the result area.

b) Anti-dengue virus IgG/IgM

- Holding the capillary tube vertically, draw the specimen (whole blood/ serum/ plasma) up to the fill line (approximately 10 µL) and add it to the specimen well (S2) of the test cassette.
- **Alternatively, a micropipette (10 µL) may be used.**
- Holding the buffer bottle vertically, add 2 drops (approximately 70 µL) of buffer to the specimen well (S2).

Avoid trapping air bubbles in the specimen well (S2) and do not add any solution to the result area.

4. Start the timer.

5. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Positive results can be visible in as little as 1 minute. Do not interpret the result after more than 20 minutes.


10. Result Interpretation

Positive for dengue virus NS1 Ag:

A coloured line develops in each control line region 'C'. Another coloured line develops in the test line region 'T'.


The test result indicates the presence of dengue virus NS1 Ag, suggesting an early, acute primary or secondary infection.

Positive for anti-dengue virus IgM:

A coloured line develops in each control line region 'C'. Another coloured line develops in the test line region 'IgM'.


The test result indicates the presence of anti-dengue virus IgM, suggesting a fresh primary dengue virus infection.

Positive for anti-dengue virus IgG:

A coloured line develops in each control line region 'C'. Another coloured line develops in the test line region 'IgG'.


The test result indicates the presence of anti-dengue virus IgG, suggesting the late stage of a primary, early secondary or previous infection.

Positive for anti-dengue virus IgM and IgG:

A coloured line develops in each control line region 'C'. A coloured line develops in the test line region 'IgM' and another one develops in the test line region 'IgG'.


The test result indicates the presence of anti-dengue virus IgM and IgG, suggesting a current primary or early secondary dengue virus infection.

Note: The colour intensity in the test line regions 'T', 'IgM' or 'IgG' may vary depending on the concentration of dengue virus NS1 antigens or anti-dengue virus antibodies in the specimen. Therefore, any shade of colour in the test line region 'T', 'IgG' or 'IgM' should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.


Negative:

A coloured line develops in each control line region 'C'. No coloured lines develop in the test line regions 'T', 'IgM' or 'IgG'.

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette:

A coloured line appearing in each control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Dengue NS1 Ag + IgG/IgM Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection and differentiation of dengue virus NS1 antigens as well as anti-dengue virus IgG and IgM in human whole blood, serum or plasma specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of dengue virus NS1 antigens or anti-dengue virus antibodies can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor may affect test results.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- The NADAL® Dengue NS1 Ag + IgG/IgM Test only detects the presence of dengue virus NS1 antigens and anti-dengue virus antibodies in specimens and should not be used as the sole criterion for a diagnosis of dengue.
- The NADAL® Dengue NS1 Ag + IgG/IgM Test cannot be used to determine dengue virus serotypes.
- Serological cross-reactivity with SARS-CoV-2 as well as across the flavivirus group (St. Louis encephalitis, West Nile virus, Japanese encephalitis and yellow fever viruses) is possible. Positive results should be confirmed using other methods. These cross-reactivities must be taken into account when making vaccination-related decisions.

- The continued presence or absence of antibodies cannot be used to determine the success or failure of therapy.
- Results from immunosuppressed patients should be interpreted with caution.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods, such as PCR or ELISA is recommended. A negative result does not at any time preclude the possibility of a recent dengue infection.
- A negative test result for dengue virus NS1 antigens can occur if the quantity of dengue virus NS1 antigens present in the specimen is below the detection limit of the test or if the detected dengue virus NS1 antigens were not present during the stage of disease in which the specimen was collected. Some patients may not produce detectable levels of IgM antibodies in early or repeated dengue infections.
- The presence of detectable dengue virus NS1 Ag may mean a positive result for an early dengue infection.
- Dengue virus NS1 antigen is expected to be detectable 1 day after the onset of fever and may persist for up to 9 days in both primary and secondary dengue infections.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® Dengue NS1 Ag + IgG/IgM Test was evaluated using clinical specimens whose status was confirmed using ELISA for the right side of the test cassette (Dengue IgG/IgM Test) and using PCR for the left side of the test cassette (Dengue NS1 Ag Test).

The results are presented in the following tables:

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgM) vs. ELISA

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgM)	Dengue IgM ELISA			Total
		Positive	Negative	
	Positive	155	5	
Negative	8	169	177	
Total	163	174	337	

Diagnostic sensitivity: 95.09% (90.62% - 97.49%)*

Diagnostic specificity: 97.13% (93.45% - 98.77%)*

Overall agreement: 96.14% (93.51% - 97.73%)*

*95% confidence interval

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG) vs. ELISA

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG)	Dengue IgG ELISA			Total
		Positive	Negative	
	Positive	138	7	
Negative	6	167	173	
Total	144	174	318	

Diagnostic sensitivity: 95.83% (91.21% - 98.08%)*

Diagnostic specificity: 95.98% (91.93% - 98.04%)*

Overall agreement: 95.91% (93.13% - 97.60%)*

*95% confidence interval

NADAL® Dengue NS1 Ag + IgG/IgM Test (NS1 Ag) vs. PCR

NADAL® Dengue NS1 Ag + IgG/IgM Test (NS1 Ag)	PCR			Total
		Positive	Negative	
	Positive	155	7	
Negative	10	192	202	
Total	165	199	364	

Diagnostic sensitivity: 93.94% (89.20% - 96.68%)*

Diagnostic specificity: 96.48% (92.92% - 98.29%)*

Overall agreement: 95.33% (92.65% - 97.06%)*

*95% confidence interval

Analytical performance

Detection limit

The detection limit of the NADAL® Dengue NS1 Ag + IgG/IgM Test is 50 ng/mL recombinant dengue virus NS1 antigen.

Interfering substances

Negative, low positive anti-dengue virus IgM and IgG serum specimens as well as low positive dengue virus NS1 Ag serum specimens spiked with the following potentially interfering substances showed no interference with the NADAL® Dengue NS1 Ag + IgG/IgM Test.

Substance	Concentration
Albumin	60 g/L
Paracetamol (Acetaminophen)	0.2 mg/mL
Atropine	0.2 mg/mL
Acetylsalicylic acid	0.2 mg/mL
Ascorbic acid	0.2 mg/mL
Bilirubin	0.2 mg/mL
Creatinine	442 µmol/L
Sodium citrat	3.80 %
Caffeine	0.2 mg/mL
EDTA	3.4 µmol/L
Haemoglobin	2 g/L
Heparin sodium	3000 IU/mL
Glucose	55 mmol/L
Salicylic acid	4.34 mmol/L

Cross-reactivity

Negative, low positive anti-dengue virus IgM and IgG serum specimens as well as low positive dengue virus NS1 Ag serum specimens were spiked with the following potentially cross-reacting substances and tested using the NADAL® Dengue NS1 Ag + IgG/IgM Test.

Substance	Concentration
anti-HAV IgM	Positive
HBsAg	Positive
anti-HCV IgM	Positive
anti- <i>H. pylori</i> IgG	Positive
anti- <i>Mycobacterium tuberculosis</i> IgM	Positive
anti- <i>Treponema pallidum</i> IgG	Positive
Rheumatoid factors (RF)	2500 IU/mL

No cross-reactivity with the specimens was observed when tested using the NADAL® Dengue NS1 Ag + IgG/IgM Test.

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive anti-dengue virus IgM and IgG serum specimens respectively, as well as the same number of negative, low, medium and high positive dengue virus NS1 Ag serum specimens using 3 lots of the NADAL® Dengue NS1 Ag + IgG/IgM tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5% - 100%). The

NADAL® Dengue NS1 Ag + IgG/IgM Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive anti-dengue virus IgM and IgG serum specimens respectively, as well as the same number of negative, low, medium and high positive dengue virus NS1 Ag serum specimens. Testing was performed by 5 operators using 3 independent NADAL® Dengue NS1 Ag + IgG/IgM test lots at 3 sites on 5 separate days. >99% of the specimens were correctly identified (375/375 correct tests per concentration, 95% confidence interval: 99.8% - 100%). The NADAL® Dengue NS1 Ag + IgG/IgM Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® Dengue NS1 Ag + IgG/IgM Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr* (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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15. Bibliographie

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German)*, 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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15. Bibliografia

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.

2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German).* 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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15. Bibliografia

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus Microbiol Immunol. 1994; 38(12): 951-5.

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15. Referências

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro*. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol*. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German)*. 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol*. 1994; 38(12): 951-5.

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15. Reference

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI. Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr* (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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