# NADAL<sup>®</sup> S. pneumoniae Test (test cassette)

REF 572004N-10

<b>DE</b>	Gebrauchsanweisung	2	<b>PL</b>	Sposób użycia	17
<b>EN</b>	Instructions for use	5	<b>PT</b>	Instruções de Utilização	20
<b>FR</b>	Instructions d'utilisation	8		Symbols	23
<b>ES</b>	Instrucciones de uso	11		Our Teams	24
<b>IT</b>	Istruzioni per l'uso	14			

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### 14. Referenzen

1. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae: virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
2. CDC. Updated recommendations for preventionof invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34) : 1102-1106, 2010.
3. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi: 10.1038/nrmicro1871, 2008.
4. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285 : 1729-1735, 2001.
5. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343 : 1917-1924, 2000.
6. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

Rev. 0, 2018-06-26 OM/KB

### 1. Intended Use

The NADAL® S. pneumoniae Test is a qualitative rapid assay which is intended to be used for the detection of *Streptococcus pneumoniae* antigen in urine without any dilution and as an aid in the diagnosis of pneumonia, meningitis and otitis media.

### 2. Introduction and Clinical Significance

*Streptococcus pneumoniae* (*S. pneumoniae*) is a gram-positive bacterium (1) first isolated from the saliva of a patient with rabies by Pasteur in 1881. The chemical structure and antigenicity of the pneumococcal capsular polysaccharide and its association with virulence and its role in human disease were explained over the period of 1915 to 1945. The bacteria are lancet-shaped anaerobic organisms spreading by direct person-to-person contact via respiratory droplets and causing serious disease in humans: The 10 most common serotypes are estimated to account for about 62% of invasive disease worldwide (2). *S. pneumoniae* colonizes upper respiratory tract tissues causing severe pneumonia and mild/ acute earache/ otitis (3). Pneumococci cause 13% to 19% of all cases of bacterial meningitis in the United States (4). One-fourth of patients with pneumococcal meningitis also have pneumonia. Clinical symptoms are generally similar to those of other forms of purulent bacterial meningitis and include headache, lethargy, vomiting, irritability, fever, nuchal rigidity, cranial signs, seizures and coma. The case-fatality of pneumococcal meningitis is about 30% but can be as high as 80% among the elderly. Bacterial pneumonia accounts for 12-16% of invasive pneumococcal disease among children aged 2 years and younger whereas *S. pneumoniae* has become the leading cause of bacterial meningitis among children younger than 5 years of age in the United States. Antibiotic treatment is efficient even if more penicillin-resistant strains have been identified (5). Several vaccines are available with variable efficiency depending on patient age or whether patients are developing some chronic illness or immunodeficiency. Nevertheless, vaccines have been demonstrated to provide protection against pneumococcal pneumonia (6).

### 3. Test Principle

The NADAL® S. pneumoniae Test is a qualitative rapid assay for the detection of *Streptococcus pneumoniae* antigen in urine without any dilution. The method employs a unique combination of monoclonal antibody-dye conjugate and monoclonal solid phase antibodies directed to different epitopes to selectively identify *Streptococcus pneumoniae* with a high degree of sensitivity. As the test sample is added, the labelled antibody-dye conjugate binds to the *S. pneumoniae* antigen forming an antibody-antigen complex. This complex binds to immobilized *S. pneumoniae* antibody in the test line region (T) producing a rose-pink-coloured line. In absence of *S. pneumoniae* in the sample, there will be no line present in the test line region (T). The reaction mixture continues flowing along the membrane.

Unbound conjugate binds to reagents in the control line region (C) producing a rose-pink-coloured line, indicating correct functioning of the test reagents and correct procedural technique.

### 4. Reagents and Materials Supplied

- 10 NADAL® S. pneumoniae test cassettes (incl. disposable pipettes)
- 1 package insert

### 5. Additional Materials Required

- Urine collection container

### 6. Storage & Stability

Test kit components should be stored in the sealed foil pouch at room temperature (4-30°C). Do not freeze the test kit. The test is stable up to the expiry date stated on the foil pouch.

### 7. Warnings and Precautions

- The test is for professional *in-vitro* diagnostic use only.
- Read the package insert carefully before performing the test.
- Do not use the test after the expiration date stated on the foil pouch.
- Do not use the test if the foil pouch is damaged.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Handle all specimens as if they contained infectious agents.
- Observe established precautions for microbiological risks throughout all procedures and standard guidelines for appropriate disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Humidity and temperature can adversely affect test results.

### 8. Specimen Collection and Preparation

Collect the urine specimen in a clean container free from detergent. If testing is not performed immediately, the specimen should be refrigerated (2-8°C) for up to 24 hours. In such case bring the specimen to room temperature prior to testing. If testing is delayed for more than 24 hours, the specimen should be frozen. The frozen specimen must be completely thawed, thoroughly mixed and brought to room temperature prior to testing. Avoid repeated freezing and thawing.

In case of doubt: It is possible to increase the sensitivity and specificity of the NADAL® S. pneumoniae Test through the following treatment of urine specimens:

- Heat a urine specimen at 100°C for 5 minutes.
  - Centrifuge it at 1000 g for 15 minutes.
  - Concentrate the supernatant using AMICON® ultra 4 centrifuge filter tubes until a 25x concentration is achieved.
- All treated/concentrated specimens should be tested following the same test procedure as for untreated specimens.

### 9. Test Procedure

**Bring samples and test cassettes to room temperature prior to testing.**

1. Remove the test cassette from the foil pouch.
2. Label the test cassette with a patient's name or control number.
3. Holding a pipette vertically, fill it with urine sample and

dispense 5 drops (200 µL) into the sample well (S) of the test cassette.


4. Read the test result after 10 minutes. Confirm negative test results after exactly 15 minutes. Do not interpret the result after more than 15 minutes.

## 10. Result Interpretation

### Positive:

In addition to the coloured line in the control line region (C) a clearly distinguishable coloured line appears in the test line region (T).

### Negative:

One coloured line appears in the control line region (C). No line appears in the test line region (T).

### Invalid:

If there is no distinct coloured line in the control line region (C), the test is invalid. It is recommended in this case that the test be repeated or a fresh urine specimen be obtained and tested later.

## 11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

## 12. Limitations

- The NADAL® S. pneumoniae Test is a screening test for the presence of *Streptococcus pneumoniae* in urine samples.
- As with any diagnostic procedure, the physician should evaluate the data obtained with the test in the light of other clinical information, including cell culture, if the results are inconsistent with the clinical presentation.
- The test does not differentiate between colonized and infected individuals.
- The NADAL® S. pneumoniae Test cannot distinguish different pneumococcal infections potentially induced by *Streptococcus pneumoniae*.
- Urine specimens with a high specific gravity (SG >1.030) may show false positive results with the NADAL® S. pneumoniae Test in rare cases.

## 13. Performance Characteristics

### Sensitivity and specificity

#### Diagnostic sensitivity and specificity

A comparative study between the BinaxNOW® Streptococcus pneumoniae rapid test and the NADAL® S. pneumoniae Test was performed using a panel of 16 fresh urine samples and 42 stored urine samples.

The NADAL® S. pneumoniae Test vs. the BinaxNOW® pneumoniae rapid test:

		BinaxNOW® Streptococcus pneumoniae rapid test (reference method)		
		+	-	Total
NADAL® S. pneumoniae Test	+	26	4	30
	-	7	21	28
Total		33	25	58

From the above table, the diagnostic sensitivity of the NADAL® S. pneumoniae Test is 79% (26/33) and the specificity is 84% (21/25) when compared to the reference method. The overall agreement is 81% (47/58) and is satisfactory despite the fact that long storage of the 42 samples (33 positive and 9 negative) might have caused a degradation of *Streptococcus pneumoniae* antigens.

### Analytical sensitivity

Serial dilutions of *Streptococcus pneumoniae* urinary antigen were performed in negative urine. The dilutions were then tested using both the NADAL® S. pneumoniae Test and BinaxNOW® test. The results presented in the table below show that the NADAL® S. pneumoniae Test is more sensitive than BinaxNOW® test: while the NADAL® S. pneumoniae Test enables the detection of the *Streptococcus pneumoniae* urinary antigen at a concentration of 35 µg/mL, the BinaxNOW® test only detects the antigen from a minimum concentration of 50 µg/mL.

Streptococcus pneumoniae urinary antigen	Results	
	NADAL® S. pneumoniae Test (after 10 minutes)	BinaxNOW® test (after 15 minutes)
25 (µg/mL)	-	-
30 (µg/mL)	-	-
35 (µg/mL)	+W	-
40 (µg/mL)	+	-
45 (µg/mL)	+	-
50 (µg/mL)	+	+W
100 (µg/mL)	+	+
200 (µg/mL)	+	+

+W = weak positive

### Cross-reactivity

Different strains of bacteria were tested using the NADAL® S. pneumoniae Test in order to determine whether they caused non-specific reactions (cross-reactions). The results are shown in the table.

Group	Result
S. pneumoniae capsular types (8)	Positive
Strep A, B, C, D, E, F, G & H	Negative
S. aureus	Negative
B. subtilis	Negative

### Intra-assay reproducibility

The intra-assay reproducibility was evaluated by running 10 replicates of 5 samples dilution having different *Streptococcus pneumoniae* concentrations. All the series showed converging results in conformance with the expected results.

Intra-lot reproducibility performance of the NADAL® S. pneumoniae Test is 100%.

#### 14. References

1. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae: virulence factors, pathogenesis, and vaccines. *Microbiological Reviews*, 59 (4), 591 – 603, 1995.
2. CDC. Updated recommendations for preventionof invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). *MMWR*, 59 (34) : 1102-1106, 2010.
3. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. *Nature Reviews, Microbiology*, 6 (4) doi: 10.1038/nrmicro1871, 2008.
4. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. *JAMA*, 285 : 1729-1735, 2001.
5. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. *N. Engl. J. Med.*, 343 : 1917-1924, 2000.
6. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. *Clin. Infect. Dis.*, 46 : 1664-1672, 2008.

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#### 14. Bibliographie

1. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J., Snippe, H. Streptococcus pneumoniae: virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
2. CDC. Updated recommendations for preventionof invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34) : 1102-1106, 2010.
3. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi: 10.1038/nrmicro1871, 2008.
4. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285 : 1729-1735, 2001.
5. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343 : 1917-1924, 2000.
6. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

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### 14. Bibliografia

1. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae: virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
2. CDC. Updated recommendations for prevention of invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34) : 1102-1106, 2010.
3. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi: 10.1038/nrmicro1871, 2008.
4. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285 : 1729-1735, 2001.
5. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343 : 1917-1924, 2000.
6. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

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

#### 14. Referências

- AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae: virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
- CDC. Updated recommendations for preventionof invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34) : 1102-1106, 2010.
- Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi: 10.1038/nrmicro1871, 2008.
- Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285 : 1729-1735, 2001.
- Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343 : 1917-1924, 2000.
- Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

Rev. 0, 2018-06-26 MHNL

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico- diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas européias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik	In-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur- begränsning	Temperatuurlimiet	Temperatur- begrænsning	Temperatur begrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tilverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

## Our Teams

### Germany:

#### Regensburg

Tel: +49 941 290 10-0  
Fax: +49 941 290 10-50

#### Moers

Tel: +49 2841 99820-0  
Fax: +49 2841 99820-1

#### Austria:

Tel: +49 941 290 10-29  
Free Tel: 0800 291 565  
Fax: +49 290 10-50  
Free Fax: 0800 298 197

#### UK & Ireland:

Tel: +49 941 290 10-18  
Free Tel –UK: 0808 234 1237  
Free Tel – IRE: 1800 555 080  
Fax: +49 290 10-50

#### France:

France Tel: 0800 915 240  
France Fax: 0800 909 493

### Switzerland

Swiss Tel: 0800 564 720  
Swiss Fax: 0800 837 476

### Belgium

Belgium Tel: 0800 718 82  
Belgium Fax: 0800 747 07

### Luxembourg

Lux. Tel: 800 211 16  
Lux. Fax: 800 261 79

### Spain:

Tel: +49 941 290 10-759  
Free Tel: 900 938 315  
Fax: +49 941 290 10-50  
Free Fax: 900 984 992

### Italy:

Tel: +49 941 290 10-34  
Fax: +49 941 290 10-50

### Poland:

Tel: +49 941 290 10-44  
Free Tel: 00 800 491 15 95  
Fax: +49 941 290 10-50  
Free Fax: 00 800 491 15 94

### Portugal:

Tel: +49 941 290 10-735  
Tel. Verde: 800 849 230  
Fax: +49 941 290 10-50  
Fax Verde: 800 849 229

### Netherlands:

Tel: +31 30 75 600  
Free Tel: 0800 0222 890  
Fax: +31 70 30 30 775  
Free Fax: 0800 024 9519

### Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607  
Free Tel: +45 80 88 87 53  
Tax: +31 703030 775

### Laboratory Diagnostics Team:

Tel: +49 941 290 10-40  
Fax: +49 941 290 10-50