# NADAL<sup>®</sup> Norovirus GI/GII+Rota- Adenovirus Test (test cassette)

REF 920007N-10

<b>DE</b> Gebrauchsanweisung	2	<b>CZ</b> Návod k použití	26
<b>EN</b> Instructions for use	6	<b>NL</b> Gebruiksaanwijzing	30
<b>FR</b> Instructions d'utilisation	10	Symbols	35
<b>ES</b> Instrucciones de uso	14	Our Teams	36
<b>IT</b> Istruzioni per l'uso	18		
<b>PL</b> Sposób użycia	22		

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### 14. Referenzen

- Parkin PC, Macarthur C, Khamballa A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. Clin. Pediatr. 49:235-239.
- Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis — An Update"
- Ciccarelli, S; Stoffi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". Infection and Drug Resistance. 6: 133–61.
- Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63–4.
- Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36–44.
- Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227–31.
- Nishio, Osamu, M. Osseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirle Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45–e80.
- Kohler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413.

Rev. 1, 2019-05-17 OM/SDe

### 1. Intended Use

The NADAL® Norovirus GI/GII+Rota-Adenovirus Test is a rapid chromatographic immunoassay for the qualitative detection and differentiation of norovirus antigens of genogroups I and II (GI and GII), as well as rotaviruses and adenoviruses in human faecal specimens. The test is intended for use as an aid in the diagnosis of norovirus, rotavirus and adenovirus infections and is designed for professional use only.

### 2. Introduction and Clinical Significance

Acute gastroenteritis is a common disorder in young children, and the associated dehydration is both a leading cause of admission to hospital in industrialised countries and a major cause of mortality in developing countries.<sup>1</sup> Symptoms may include diarrhoea, vomiting and abdominal pain.<sup>2</sup> Fever, a lack of energy and dehydration may also occur.<sup>3</sup> This typically lasts less than two weeks.<sup>4</sup> Enteric viruses have been recognised as the most significant etiological agents of the disease, with rotavirus, adenovirus 40/41, norovirus and astrovirus currently being recognised as the most significant pathogens of childhood viral gastroenteritis.<sup>5,6</sup>

Rotavirus is the most common cause of gastroenteritis in children.<sup>7</sup> Viruses cause approximately 70% of episodes of infectious diarrhoea in the paediatric age group.<sup>8</sup>

Adenoviruses, primarily Ad40 and Ad41, are a leading cause of diarrhoea in children, second only to rotaviruses. Infections are most frequently observed in children under two years old, but have been found in patients of all ages. Further studies indicate that adenoviruses are associated with 4-15% of all hospitalised cases of viral gastroenteritis.<sup>9</sup>

Norovirus is a leading cause of gastroenteritis among adults in United States of America, causing more than 90% of outbreaks.<sup>10</sup> These localised epidemics typically occur when groups of people spend time in close physical proximity to each other, such as on cruise ships, in hospitals or in restaurants. Noroviruses commonly isolated in cases of acute gastroenteritis belong to two genogroups: GI and GII. Patients may remain infectious even after their diarrhoea has cleared up.<sup>10</sup>

Gastroenteritis is typically diagnosed clinically, based on a person's signs and symptoms.<sup>10</sup> However, stool cultures should be performed in those with blood in their stool, those who may have been exposed to food poisoning and in those who have recently traveled to developing countries.<sup>8</sup> They may also be appropriate for children under 5 years old, elderly people, and those with poor immune function.<sup>11</sup> Diagnostic testing may also be carried out for monitoring.<sup>10</sup>

### 3. Test Principle

The NADAL® Norovirus GI/GII+Rota-Adenovirus Test enables the qualitative detection of norovirus antigens GI and GII, as well as rotaviruses and adenoviruses, through the visual interpretation of colour development on the internal test strips.

The test cassette contains three different test strips for the qualitative detection of norovirus antigens GI/GII, rotaviruses and adenoviruses.

Anti-norovirus GI/GII, anti-rotavirus and anti-adenovirus antibodies are immobilised in the test line regions (GI) and

(GII) for noroviruses GI and GII, as well as in the test line regions (T) for rotaviruses and adenoviruses of the membrane. During the test, extracted antigens bind to anti-norovirus GI/GII, anti-rotavirus and anti-adenovirus antibodies, which are conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-norovirus GI/GII, anti-rotavirus and anti-adenovirus antibodies in the test line regions (GI), (GII) and (T) respectively. Excess coloured particles are captured in the control line regions (C).

The presence of a red line in the test line region (GI) and/or (GII) and/or (T) indicates a positive result for particular viral antigens, whereas its absence indicates a negative result. The formation of a blue line in each control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

### 4. Reagents and Materials Supplied

- 10 NADAL® Norovirus GI/GII+Rota-Adenovirus test cassettes (incl. disposable pipettes for liquid faecal samples)
- 10 specimen collection tubes with buffer
- 1 package insert

### 5. Additional Materials Required

- Specimen collection containers
- Disposable gloves
- Timer

### 6. Storage & Stability

NADAL® Norovirus GI/GII+Rota-Adenovirus test kits should be stored at 2-30°C and used by the expiry date indicated on the packaging. Test cassettes should remain in sealed pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date. Care should be taken to protect components of the test kit from contamination. Do not use the test if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Failure to bring specimens and reagents to room temperature prior to testing may decrease assay sensitivity.

Inaccurate or inappropriate specimen collection, storage or transport may yield false negative test results.

- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures as well as standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

## 8. Specimen Collection and Preparation

### Specimen Collection and Storage:

Use clean, dry specimen collection containers for specimen collection. The best results will be obtained if the assay is performed within 2 hours of specimen collection.

**Note:** Specimens collected in the specimen collection container may be stored at 2-8°C for up to 1-2 days, or at -20°C for up to 3 months if they are not to be tested within 2 hours of preparation.

Viral detection can be improved by collecting specimens at the onset of symptoms. If specimens are collected long after the onset of diarrhoeic symptoms, either the antigen quantity may not be sufficient enough to obtain a positive reaction or the antigens detected may not be linked to the diarrhoeic episode.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

### Specimen preparation:

Use a separate specimen collection tube with buffer for each sample.


Hold the specimen collection tube upright and remove the applicator stick by unscrewing the light-blue cap. Be careful not to spill or splash buffer solution from the collection tube.

### For solid specimens:

Collect approx. 50 mg of faeces (equivalent to ¼ of a pea) by inserting the applicator stick in at least 3 different sites of the faecal specimen.

### For liquid specimens:

Holding the pipette vertically, collect a faecal specimen and add 2 drops (approx. 80 µL) to the specimen collection tube containing buffer.

Place the applicator stick, along with the faecal specimen (or without the faecal specimen if the specimen is liquid), back into the specimen collection tube and tightly screw the cap closed.

Shake the specimen collection tube to thoroughly mix the specimen and buffer.

Specimens prepared in the specimen collection tubes can be stored at -20°C for up to 6 months if not tested within 1 hour after preparation.

## 9. Test Procedure

**Bring tests and samples with buffer to room temperature (15-30°C) prior to testing. Do not open the foil pouch until you are ready to perform the assay.**


1. Remove the NADAL® Norovirus GI/GII+Rota-Adenovirus test cassette from the sealed foil pouch and use it as soon as possible. Label the test cassette with the patient or control identification. For best results, the assay should be performed within 1 hour.
2. Shake the specimen collection tube to ensure good sample dispersion. Unscrew the white cap and break off the tip of the tube using a piece of tissue. Holding the specimen collection tube vertically, dispense 2 drops of solution into each sample well (S) of the test cassette.

**Avoid trapping air bubbles in the sample wells (S) and do not add any solution to the result areas.**

3. Start the timer.

**Note:** If the specimen does not migrate due to the presence of particles, centrifuge the specimen contained in the specimen collection tube. Collect 80 µL of supernatant, dispense it into each specimen well (S) of a new test cassette and, following the instructions stated above, start the test again.


4. Wait for the coloured line(s) to appear. Read the test result after 10 minutes. Do not interpret the result after more than 20 minutes.

## 10. Result Interpretation


### Positive for norovirus GI (Noro)

One blue line develops in the control line region (C) and one red line develops in the test line region (GI) for norovirus GI (Noro).

### Positive for norovirus GII (Noro)


One blue line develops in the control line region (C) and one red line develops in the test line region (GII) for norovirus GII (Noro).


**Positive for rotavirus (Rota):**

One blue line develops in the control line region (C) and one red line develops in the test line region (T) for rotavirus (Rota).

**Positive for adenovirus (Adeno):**

One blue line develops in the control line region (C) and one red line develops in the test line region (T) for adenovirus (Adeno).


**Note:** A combination of positive results for various parameters is possible.


**Note:** The colour intensity in the test line regions may vary depending on the concentration of the analyte present in the specimen. Therefore, any shade of colour in the test line regions should be considered positive.

**Negative:**

One blue line develops in each control line region (C). No lines develop in test line regions (GI and (GII) and in test line regions (T).

**Invalid:**

A blue control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

**11. Quality Control**

An internal procedural control is included in the test cassette:

A blue line appearing in each control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP)* recommends the use of external control materials to ensure proper test kit performance.

**12. Limitations**

- The NADAL® Norovirus GI/GII+Rota-Adenovirus Test is for professional *in-vitro* diagnostic use only and should only be used for the qualitative detection of noroviruses GI and GII, as well as rotaviruses and adenoviruses. The colour intensity of the test lines cannot be used for semi-quantitative or quantitative evaluation.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

- The sections 'Test Procedure' and 'Result Interpretation' must be followed closely while testing. Failure to follow the procedure may lead to inaccurate test results.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a norovirus, rotavirus and/or adenovirus infection, as the concentration of viral particles may be below the detection limit of the test.
- In order to prevent pseudo-outbreaks, positive test results in newborns must be confirmed using an alternative testing method (PCR).
- If visible blood is present in faecal specimens, false-positive results cannot be ruled out and must therefore be verified using an alternative test method.

**13. Performance Characteristics****Clinical Performance**

A total of 210 faecal samples from individual subjects were tested using the NADAL® Norovirus GI/GII+Rota-Adenovirus Test (adenovirus) and a commercially available ELISA.

The results are presented in the following table:

**NADAL® Norovirus GI/GII+Rota-Adenovirus Test (adenovirus) vs. ELISA**

		NADAL® Norovirus GI/GII+Rota-Adenovirus (adenovirus)		
		+	-	Total
ELISA	+	82	1	83
	-	0	127	127
	Total	82	128	210

Relative sensitivity: 98.8% (93.5% - 99.8%)\*

Relative specificity: >99.9% (97.1% - 100.0%)\*

Overall agreement: 99.5% (97.4% - 99.9%)\*

\*95% Confidence interval

A total of 242 faecal samples from individual subjects were tested using the NADAL® Norovirus GI/GII+Rota-Adenovirus Test (rotavirus) and a commercially available ELISA.

The results are presented in the following table.

**NADAL® Norovirus GI/GII+Rota-Adenovirus Test (rotavirus) vs. ELISA:**

		NADAL® Norovirus GI/GII+Rota-Adenovirus Test (rotavirus)		
		+	-	Total
ELISA	+	79	3	82
	-	0	160	160
	Total	79	163	242

Relative sensitivity: 96.3% (89.8% - 98.7%)\*

Relative specificity: >99.9% (97.7% - 100.0%)\*

Overall agreement: 98.8% (96.4% - 99.6%)\*

\*95% Confidence interval

A total of 398 faecal samples were tested using the NADAL® Norovirus GI/GII+Rota-Adenovirus Test and a commercially available immunochromatographic test for the detection of norovirus antigens of genogroups I and II.

The results are presented in the following table.

**NADAL® Norovirus GI/GII+Rota-Adenovirus Test (norovirus) vs. another Norovirus GI/GII rapid test:**

14. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. *American Journal of Infection Control* 2018; 46(12):1411–1413.

Rev. 1, 2019-05-17 OM/SDe

		NADAL® Norovirus GI/GII+Rota-Adenovirus Test (norovirus)		
		+	–	Total
Another Norovirus GI/GII rapid test	+	142	3	145
	–	1	252	253
	Total	143	255	398

Relative sensitivity: 97.9% (94.1% - 99.3%)\*

Relative specificity: 99.6% (97.8% - 99.9%)\*

Overall agreement: 99.0% (97.4% - 99.6%)\*

\*95% Confidence interval

**Cross-reactivity:**

Cross-reactivity with the following organisms was studied at a concentration of  $1.0 \times 10^9$  organisms/mL. Faecal specimens containing the following organisms were found to be negative when tested using the NADAL® Norovirus GI/GII+Rota-Adenovirus Test.

<i>Staphylococcus aureus</i>	<i>Enterococcus faecalis</i>
<i>Acinetobacter spp.</i>	<i>Moraxella catarrhalis</i>
<i>Proteus vulgaris</i>	<i>Chlamydia trachomatis</i>
<i>Enterococcus faecium</i>	<i>Pseudomonas aeruginosa</i>
<i>Gardnerella vaginalis</i>	<i>Neisseria gonorrhoeae</i>
<i>Klebsiella pneumoniae</i>	<i>Group C Streptococcus</i>
<i>Acinetobacter calcoaceticus</i>	<i>Helicobacter pylori</i>
<i>Proteus mirabilis</i>	<i>Salmonella Choleraesuis</i>
<i>Group B Streptococcus</i>	<i>Haemophilus influenzae</i>
<i>Escherichia coli</i>	<i>Neisseria meningitidis</i>
<i>Candida albicans</i>	

**14. References**

- Parkin PC, Macarthur C, Khambalia A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. *Clin. Pediatr.* 49:235-239.
- Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis – An Update"
- Ciccarelli, S; Stolfi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". *Infection and Drug Resistance.* 6: 133–61.
- Schlossberg, David (2015). *Clinical infectious disease* (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". *Recent Patents on Anti-infective Drug Discovery.* 6 (1): 54–63.
- Dennehy PH (January 2011). "Viral gastroenteritis in children". *The Pediatric Infectious Disease Journal.* 30 (1): 63–4.
- Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". *Current Opinion in Gastroenterology.* 26 (1): 36–44.
- Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". *Australian Family Physician.* 34 (4): 227–31.
- Nishio, Osamu, M. Ooseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." *Microbiol. Immunol.* 1990; 34(10): 871-877.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". *Recent Patents on Anti-infective Drug Discovery.* 6 (1): 54–63.
- Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirle Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". *Clinical Infectious Diseases.* 65: e45–e80.
- Kohler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. *J Pediatr Gastroenterol Nutr* 2008; 46(4):471-472.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. *Journal of Microbiology, Immunology and Infection* 2016; 49:947e954.

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**14. Bibliographie**

- Parkin PC, Macarthur C, Khambalia A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. Clin. Pediatr. 49:235-239.
- Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis — An Update"
- Ciccarelli, S; Stolfi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". Infection and Drug Resistance. 6: 133-61.
- Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54-63.
- Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63-4.
- Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36-44.
- Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227-31.
- Nishio, Osamu, M. Osseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54-63.
- Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirlé Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45-e80.
- Köhler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413.

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**14. Referencias**

1. Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
2. Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
3. Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63–4.
4. Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36–44.
5. Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227–31.
6. Nishio, Osamu, M. Oseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
7. Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
8. Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirle Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45–e80.
9. Köhler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
10. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
11. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413.

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#### 14. Bibliografia

- Parkin PC, MacArthur C, Khambalia A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. Clin. Pediatr. 49:235-239.
- Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis — An Update"
- Ciccarelli, S; Stolfi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". Infection and Drug Resistance. 6: 133–61.
- Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63–4.
- Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36–44.
- Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227–31.
- Nishio, Osamu, M. Ooseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirlé Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45–e80.
- Köhler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413.

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#### 14. Bibliografia

1. Parkin PC, Macarthur C, Khambalia A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. Clin. Pediatr. 49:235-239.
2. Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis — An Update"
3. Ciccarelli, S; Stolfi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". Infection and Drug Resistance. 6: 133-61.
4. Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
5. Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54-63.
6. Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63-4.
7. Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36-44.
8. Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227-31.
9. Nishio, Osamu, M. Osseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
10. Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54-63.
11. Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirlé Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45-e80.
12. Kohler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
13. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
14. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413.

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#### 14. Reference

- Parkin PC, Macarthur C, Khambalia A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. Clin. Pediatr. 49:235-239.
- Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis – An Update"
- Ciccarelli, S; Stolfi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". Infection and Drug Resistance. 6: 133–61.
- Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63–4.
- Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36–44.
- Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227–31.
- Nishio, Osamu, M. Ooseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirle Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45–e80.
- Kohler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413.

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



**14. Referenties**

- Parkin PC, Macarthur C, Khambalia A, Goldman RD, Friedman JN. June 2009. Clinical and laboratory assessment of dehydration severity in children with acute gastroenteritis. Clin. Pediatr. 49:235-239.
- Singh, Amandeep (July 2010). "Pediatric Emergency Medicine Practice Acute Gastroenteritis — An Update"
- Ciccarelli S, Stolfi, I; Caramia, G (29 October 2013). "Management strategies in the treatment of neonatal and pediatric gastroenteritis". Infection and Drug Resistance. 6: 133–61.
- Schlossberg, David (2015). Clinical infectious disease (Second ed.). p. 334. ISBN 9781107038912. Archived from the original on 2017-09-08.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Dennehy PH (January 2011). "Viral gastroenteritis in children". The Pediatric Infectious Disease Journal. 30 (1): 63–4.
- Szajewska, H; Dziechciarz, P (January 2010). "Gastrointestinal infections in the pediatric population". Current Opinion in Gastroenterology. 26 (1): 36–44.
- Webb, A; Starr, M (April 2005). "Acute gastroenteritis in children". Australian Family Physician. 34 (4): 227–31.
- Nishio, Osamu, M. Osseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. "Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces." Microbiol. Immunol. 1990; 34(10): 871-877.
- Eckardt AJ, Baumgart DC (January 2011). "Viral gastroenteritis in adults". Recent Patents on Anti-infective Drug Discovery. 6 (1): 54–63.
- Shane, Andi L; Mody, Rajal K; Crump, John A; Tarr, Phillip I; Steiner, Theodore S; Kotloff, Karen; Langley, Joanne M; Wanke, Christine; Warren, Cirle Alcantara; Cheng, Allen C; Cantey, Joseph; Pickering, Larry K (19 October 2017). "2017 Infectious Diseases Society of America Clinical Practice Guidelines for the Diagnosis and Management of Infectious Diarrhea". Clinical Infectious Diseases. 65: e45–e80.
- Kohler H, Jungert J, Korn K. Norovirus pseudo-outbreak in a neonatal intensive care unit. J Pediatr Gastroenterol Nutr 2008; 46(4):471-472.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954.
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413.

Rev. 1, 2019-05-17 HB

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-mærking	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik	In-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur-begränsning	Temperatuurlimiet	Temperatur-begrænsning	Temperatur begrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertäkäyttöinen	Får inte återvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il l ingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrekkelig for <n> tester

## Our Teams

### Germany:

#### Regensburg

Tel: +49 941 290 10-0  
 Fax: +49 941 290 10-50

#### Moers

Tel: +49 2841 99820-0  
 Fax: +49 2841 99820-1

#### Austria:

Tel: +49 941 290 10-29  
 Free Tel: 0800 291 565  
 Fax: +49 290 10-50  
 Free Fax: 0800 298 197

#### UK & Ireland:

Tel: +49 941 290 10-18  
 Free Tel – UK: 0808 234 1237  
 Free Tel – IRE: 1800 555 080  
 Fax: +49 290 10-50

#### France:

France Tel: 0800 915 240  
 France Fax: 0800 909 493

### Switzerland

Swiss Tel: 0800 564 720  
 Swiss Fax: 0800 837 476

### Belgium

Belgium Tel: 0800 718 82  
 Belgium Fax: 0800 747 07

### Luxembourg

Lux. Tel: 800 211 16  
 Lux. Fax: 800 261 79

### Spain:

Tel: +49 941 290 10-759  
 Free Tel: 900 938 315  
 Fax: +49 941 290 10-50  
 Free Fax: 900 984 992

### Italy:

Tel: +49 941 290 10-34  
 Fax: +49 941 290 10-50

### Poland:

Tel: +49 941 290 10-44  
 Free Tel: 00 800 491 15 95  
 Fax: +49 941 290 10-50  
 Free Fax: 00 800 491 15 94

### Portugal:

Tel: +49 941 290 10-735  
 Tel. Verde: 800 849 230  
 Fax: +49 941 290 10-50  
 Fax Verde: 800 849 229

### Netherlands:

Tel: +31 30 75 600  
 Free Tel: 0800 0222 890  
 Fax: +31 70 30 30 775  
 Free Fax: 0800 024 9519

### Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607  
 Free Tel: +45 80 88 87 53  
 Tax: +31 703030 775

### Laboratory Diagnostics Team:

Tel: +49 941 290 10-40  
 Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH