# NADAL® Strep B Test (test cassette)

REF 232001

|           |                            |    |           |                 |    |
|-----------|----------------------------|----|-----------|-----------------|----|
| <b>DE</b> | Gebrauchsanweisung         | 2  | <b>CZ</b> | Návod k použití | 31 |
| <b>EN</b> | Instructions for use       | 7  | <b>FI</b> | Käyttöohje      | 35 |
| <b>FR</b> | Instructions d'utilisation | 11 |           | Symbols         | 39 |
| <b>ES</b> | Instrucciones de uso       | 15 |           | Our Teams       | 40 |
| <b>IT</b> | Istruzioni per l'uso       | 19 |           |                 |    |
| <b>PL</b> | Sposób użycia              | 23 |           |                 |    |
| <b>PT</b> | Instruções de Utilização   | 27 |           |                 |    |

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### 15. Referenzen

- Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976
- Lim, D.V., Morales, W.J., Walsh, W.J. and Kazanis, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986
- Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910
- Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336-

5. Hahn H., Kaufmann S.H.E., Schulz T.F., Suerbaum S. (eds) Medizinische Mikrobiologie, Springer Verlag, 6. komplett überarbeitete Auflage, Chapter 26.2 Streptococcus agalactiae
6. Verani J.R., McGee L., Schrag S.J. (2010) Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010

Rev. 1, 2017-05-30 OM/UJ

### 1. Intended Use

The NADAL® Strep B Test is a rapid, visual immunoassay for the qualitative, presumptive detection of group B Streptococcus (GBS) antigens in specimens collected from vaginal or rectal swab. Ear or throat swabs from newborns collected directly after birth can also be used. The test is intended as an aid in diagnosing the presence of GBS in sample material that may either be a sign of colonisation or infection. This test provides epidemiological information about infections caused by GBS. The NADAL® Strep B Test is designed for professional use only. The detection limit of the assay is  $1 \times 10^5$  bacteria/swab and does not depend on viability of organisms.

### 2. Introduction and Clinical Significance

Group B Streptococci (GBS) or *Streptococcus agalactiae* are among the most frequent causes of life-threatening infections in neonates. Up to 30% of all pregnant women are colonised with GBS. Some recent studies have shown that the intrapartum antibiotic prophylaxis of GBS-colonised women significantly reduces the incidence of GBS-caused sepsis in newborns. Routine examination for GBS is frequently recommended between the 35<sup>th</sup> and the 37<sup>th</sup> week of pregnancy. A CDC study has shown that this screening approach is 50% more effective than the use of antibiotics for pregnant women identified by clinical risk approach.

Standard culture methods require 24 to 48 hours, and the results may not be available soon enough for efficient treatment. Thus, methods utilising more rapid screening techniques are required. The NADAL® Strep B Test is especially suitable if time and/or availability of culture methods are limited. It enables the detection of GBS directly from swabs, helping physicians make a presumptive diagnosis and decide whether a therapy should be indicated.

### 3. Test Principle

The NADAL® Strep B Test detects group B Streptococcus antigens through visual interpretation of colour development on the internal test strip. Anti-streptococcus B antibodies are immobilised in the test line region of the membrane. During testing, the specimen reacts with further polyclonal anti-streptococcus B antibodies conjugated to coloured particles and precoated onto the conjugate pad of the internal test strip. The mixture then migrates along the membrane by capillary action and interacts with reagents on the membrane.

If there are sufficient streptococcus B antigens in the specimen, a coloured line will form in the test line region of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The appearance of a coloured line in the control line region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

### 4. Reagents and Materials Supplied

- 20 NADAL® Strep B test cassettes (individually packaged). Each test contains coloured conjugates and reactive reagents precoated in the corresponding regions of the internal test strip.
- 1 Reagent 1: 1.0 M sodium nitrite

Danger

H301: Toxic if swallowed

- 1 Reagent 2: 0.4 M acetic acid
- 1 Positive Control: Non-viable GBS: 0,09% sodium azide
- 1 Negative Control: containing heat-killed non-Group B Streptococci
- Provided additional material according to 93/42/EEC:  
20 sterilized swabs Copan Floqswabs: For specimen collection; CE 0123

Copan Italia S.p.A.,  
Via Perotti 10, 25125 Brescia, Italy

- 20 disposable pipettes: For adding specimens
- 20 extraction tubes: For specimen preparation
- 1 reagent holder
- 1 package insert

### 5. Additional Materials Required

- Timer

### 6. Storage & Stability

The test kit should be stored at 2-30°C and used before the expiry date printed on the packaging. The test cassette should remain in the sealed foil pouch until use. Do not freeze test kits. Care should be taken to protect components of the test kit from contamination. Do not use tests if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not swap caps between different extraction reagent bottles.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled

in accordance with usual safety precautions (e.g., do not ingest or inhale).

- Use the provided swabs for sample collection. Do not use calcium alginate, cotton-tipped or wooden-shafted swabs.
- Do not use swabs from damaged pouches.
- Reagents 1 & 2 are slightly caustic. Avoid contact with eyes or mucous membranes. In the event of accidental contact, wash thoroughly with water.
- Controls contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of this solution, always flush with copious amount of water to prevent azide build-up. Avoid contact with eyes or mucous membranes. In the event of accidental contact, wash thoroughly with water.
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

## 8. Specimen Collection and Preparation

- The quality of specimen obtained is of extreme importance. Collect swab specimens using standard clinical methods.
- It is recommended that swab specimens be processed as soon as possible after collection. If swabs are not processed immediately, they should be placed into a sterile, dry, tightly capped tube or bottle and refrigerated. Do not freeze swabs. Swabs can be stored at room temperature (15-30°C) up to 4 hours or refrigerated (2-8°C) up to 24 hours. All specimens should be brought to room temperature (15-30°C) before testing.
- Do not place swabs in any transport device containing liquid transport media or transport media containing agar or charcoal. Transport media may interfere with the assay and viability of organisms. If transport medium is required, we recommend using Modified Stuart's Transport Medium as outlined in the manufacture's instructions.
- If a bacterial culture is desired, lightly roll the swab on an appropriate cell culture plate before using it in the test. The extraction reagents in the test will kill bacteria on swabs and make them impossible to culture.

## 9. Test Procedure

**Bring tests, specimens and reagents to room temperature (15-30°C) prior to testing.**

**To avoid cross contamination, do not allow the tips of the reagent bottles to come into contact with sample material.**

**For testing of controls refer to section 11 'Quality Control'.**

### 1. Prepare swab specimens:

- Place a clean extraction tube onto the designated area of the reagent holder. Add 4 drops of Reagent 1 and 4 drops of Reagent 2 to the extraction tube. In order to ensure reliable drop size when adding the reagents, hold the dropper bottles vertically. Mix the solution by gently swirling the extraction tube.
- Immediately immerse the swab with the sample into the extraction tube. Using circular motions, roll the swab against the side of the extraction tube and squeeze the extraction tube tightly so that as much liquid

as possible is expressed from the swab and can reabsorb. Repeat at least 5 times.

- Let the solution stand for 3-5 minutes at room temperature. Then express the liquid from the swab head by rolling the swab against the wall of the extraction tube and squeezing the extraction tube as the swab is withdrawn. Discard the swab in accordance with the guidelines for handling infectious agents.
- 2. Remove the test cassette from the sealed foil pouch and place it on a clean and level surface. The test cassette should be used immediately or within one hour at the latest after opening the foil pouch. Label the test cassette with the patient or control identification.
- 3. Transfer 3 drops (approximately 120 µL) of the extracted solution or the control with the included disposable pipette from the extraction tube to the sample well of the test cassette.

**Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.**

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

- 4. Wait for the coloured line(s) to appear. The test result should be read after 10 minutes. Do not interpret the result after more than 20 minutes.

## 10. Result Interpretation

### Positive

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T). This indicates that streptococcus B antigen has been detected in the sample.

### Negative

Only one coloured line appears in the control line region (C). No apparent coloured line appears in the test line region (T). This indicates that no streptococcus B antigen has been detected or the GBS concentration is below the detection limit.

### Invalid

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

**Note:**

The colour intensity in the test line region (T) may vary depending on the concentration of the antigen present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of the antigen in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

**11. Quality Control**

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP)* recommends the use of control materials to ensure proper test kit performance. A positive control containing heat-killed group B streptococci and a negative control containing heat-killed non-group B streptococci are provided with each test kit.

**Operating Procedure for External Quality Control Testing:**

**Bring tests, reagents and controls to room temperature (15-30°C) prior to testing.**

1. Place a clean extraction tube onto the designated area of the reagent holder. Add 4 drops of Reagent 1 and 4 drops of Reagent 2 to the extraction tube. In order to ensure reliable drop size when adding the reagents, hold the dropper bottles vertically. Mix the solution by gently swirling the extraction tube.
2. Thoroughly mix the control solution (Positive/Negative Control) by shaking the bottle vigorously. Add 1 drop of Positive or Negative Control to the tube.
3. Place a clean sterile swab into the tube and swirl. Leave the swab in the extraction tube for 1 minute. Then express the liquid from the swab head by rolling the swab against the wall of the extraction tube and squeezing the extraction tube as the swab is withdrawn. Discard the swab.
4. Continue as described in the Step 2 of the 'Test Procedure'. If controls do not yield expected results, do not use the tests with samples. Repeat the Quality Control Testing or contact your distributor.

**12. Limitations**

- The NADAL® Strep B Test is for professional *in-vitro* diagnostic use and should only be used for the qualitative detection of group B streptococci (GBS). No meaning should be inferred from the colour intensity or width of any apparent lines.
- The accuracy of the test depends on the quality of the swab specimen. False negative results may occur due to improper specimen collection or storage. A negative result may also be obtained from patients with a light GBS colonization due to low antigen concentration. It is advisable to confirm negative test results with an alternative method e.g. culture.
- Maternal colonisation with GBS can be intermittent, transient or permanent and usually does not cause any

clinical symptoms. Therefore, it is usually recommended that the time point of testing should be close to birth.

- For throat/ear swabs testing of newborns no clinical studies are available for the test. Negative results should not be used to exclude the presence of GBS because bacterial numbers may be below the detection limit. It is generally recommended to take swabs from multiple locations to increase the probability of detection.
- The test does not differentiate asymptomatic carriers of Group B streptococci from those with infection. If clinical signs and symptoms are not consistent with laboratory test results, a follow-up culture is recommended.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

**13. Expected values**

Up to 30% of women carry GBS in their vaginal or rectal flora. Usually this colonisation does not cause any symptoms. However, severe infections may occur if bacteria are transmitted during birth so that GBS is a leading cause of morbidity and mortality amongst newborns. The incidence of early-onset neonatal disease is approximately one to two cases per 1,000 live births. Mortality rate is up to 20% in affected neonates. In heavily colonized women infection risk is increased ca. 2.5 fold if compared to lightly colonised women. Intrapartum chemoprophylaxis reduces the infection risk of early onset neonatal disease.

**14. Performance Characteristics****Correlation Study**

A correlation study between the NADAL® Strep B Test and conventional culture was performed. Vaginal, rectal swabs were used as sample material.

The results are presented in the following table:

| NADAL®<br>Strep B<br>Test |       | Culture |    |       |
|---------------------------|-------|---------|----|-------|
|                           |       | +       | –  | Total |
|                           | +     | 20      | 2  | 22    |
|                           | –     | 2       | 96 | 98    |
|                           | Total | 22      | 98 | 120   |

Relative sensitivity: 90.9%

Relative specificity: 97.9%

Overall agreement: 96.7%

**Detection limit**

The detection limit of the NADAL® Strep B Test was determined to be  $1 \times 10^5$  organisms/swab. Sensitivity testing is done for each lot as part of quality control and is stated in the lot-specific Certificate of Analysis.

**Prozone effect study**

No adverse effect on T-line formation was recorded for streptococcus B concentration of up to  $1 \times 10^9$  organisms/swab.

**Specificity study**

Cross-reactivity studies with organisms likely to be found in the vagina/rectum or with other *Streptococcus* species were performed using the NADAL® Strep B Test. The following

organisms were tested at  $1 \times 10^7$  organisms/swab and showed negative results:

6. Verani J.R., McGee L., Schrag S.J. (2010) Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010

Rev. 1, 2017-05-30 OM/UJ

| Organism                           | ATCC No. | Organism                         | ATCC No. |
|------------------------------------|----------|----------------------------------|----------|
| <i>Candida albicans</i>            | 1106     | Strep C                          | 12401    |
| <i>Corynebacterium diphtheriae</i> | 13812    | Strep F                          | 12392    |
| <i>Enterococcus durans</i>         | 19432    | Strep G                          | 12394    |
| <i>Enterococcus faecalis</i>       | 19433    | <i>Streptococcus canis</i>       | 43496    |
| <i>Haemophilus influenzae</i>      | 9006     | <i>Streptococcus equisimilis</i> | 9528     |
| <i>Klebsiella pneumoniae</i>       | 9987     | <i>Streptococcus equisimilis</i> | 9542     |
| <i>Neisseria gonorrhoeae</i>       | 27633    | <i>Streptococcus equisimilis</i> | 12388    |
| <i>Neisseria meningitidis</i>      | 13077    | <i>Streptococcus mutans</i>      | 25175    |
| <i>Neisseria sicca</i>             | 9913     | <i>Streptococcus pneumoniae</i>  | 27338    |
| <i>Neisseria subflava</i>          | 14799    | <i>Streptococcus sanguis</i>     | 10556    |
| <i>Pseudomonas aeruginosa</i>      | 9721     | <i>Streptococcus oralis</i>      | 9811     |
| <i>Serratia marcescens</i>         | 8100     | <i>Streptococcus mitis</i>       | 903      |
| <i>Staphylococcus aureus</i>       | 12598    | <i>Streptococcus anginosus</i>   | 33397    |
| <i>Staphylococcus epidermidis</i>  | 1228     | <i>Streptococcus intermedius</i> | 27335    |
| Strep A                            | 19615    |                                  |          |

#### Physician Office Laboratory (POL) Studies

An evaluation of the NADAL® Strep B Test was conducted at three physicians' office laboratory sites, using a panel of coded samples containing negative control, low positive and medium positive specimens. Each specimen level was tested at each site in replicates of 30 over a period of 3 days. The study showed >99.9% agreement with the expected results.

#### Interference study

A variety of vaginal washes were tested at concentrations of 1%. None of them interfered with the generation of correct test results.

#### Inter-lot and intra-lot variability

Three independent lots were tested with negative, low, medium and high positive controls in 10-fold determinations. No unexpected or inconsistent results were obtained, indicating that inter-lot and intra-lot variability is low.

#### 15. References

1. Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976
2. Lim, D.V., Morales, W.J., Walsh, W.J. and Kazanis, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986
3. Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910
4. Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336-
5. Hahn H., Kaufmann S.H.E., Schulz T.F., Suerbaum S. (eds) Medizinische Mikrobiologie, Springer Verlag, 6. komplett überarbeitete Auflage, Chapter 26.2 Streptococcus agalactiae

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### 15. Bibliographie

- Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976
- Lim, D.-V., Morales, W.J., Walsh, W.J. and Kazanin, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986
- Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910
- Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336-
- Hahn H., Kaufmann S.H.E., Schulz T.F., Suerbaum S. (eds) Medizinische Mikrobiologie, Springer Verlag, 6. komplett überarbeitete Auflage, Chapter 26.2 Streptococcus agalactiae
- Verani J.R., McGee L., Schrag S.J. (2010) Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010

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15. Referencias

1. Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976
2. Lim, D.V., Morales, W.J., Walsh, W.J. and Kazanis, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986
3. Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910
4. Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336
5. Hahn H., Kaufmann S.H.E., Schulz T.F., Suerbaum S. (eds) Medizinische Mikrobiologie, Springer Verlag, 6. komplett überarbeitete Auflage, Chapter 26.2 Streptococcus agalactiae
6. Verani J.R., McGee L., Schrag S.J. (2010) Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010

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### 15. Bibliografia

- Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976
- Lim, D.V., Morales, W.J., Walsh, W.J. and Kazanis, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986
- Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910
- Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336-
- Hahn H., Kaufmann S.H.E., Schulz T.F., Suerbaum S. (eds) Medizinische Mikrobiologie, Springer Verlag, 6. komplett überarbeitete Auflage, Chapter 26.2 Streptococcus agalactiae
- Verani J.R., McGee L., Schrag S.J. (2010) Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010

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#### 15. Bibliografia

1. Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976
2. Lim, D.V., Morales, W.J., Walsh, W.J. and Kazanis, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986
3. Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910
4. Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336-
5. Hahn H., Kaufmann S.H.E., Schulz T.F., Suerbaum S. (eds) Medizinische Mikrobiologie, Springer Verlag, 6. komplett überarbeitete Auflage, Chapter 26.2 Streptococcus agalactiae
6. Verani J.R., McGee L., Schrag S.J. (2010) Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010

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15. Referências

1. Finch, R.G., French, G.L., and Phillips, I.; Group B streptococci in the female genital tract; Br. Med. J., 1 (6020) 1245-1247, 1976

2. Lim, D.V., Morales, W.J., Walsh, W.J. and Kazanin, D.; Reduction of morbidity and mortality rates for neonatal group B streptococcal disease through early diagnosis and chemo-prophylaxis; J. Clin. Microbiol. 23 489-492, 1986

3. Apgar B.S. et al. (2005) Prevention of Group B Streptococcal Disease in the Newborn, American Family Physician 71(5): 903-910

4. Ippolito D.L. et al. (2010) Group B Streptococcus serotype prevalence in reproductive-age women at a tertiary care military medical center relative to global serotype distribution, BMC Infectious Diseases 10: 336-

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nal von minden GmbH