NADAL[®] Influenza A/B Test (test cassette)

REF 242015N-10

DE Gebrauchsanweisung	2	CZ Návod k použití	26
EN Instructions for use	6	FI Käyttöohje	30
FR Instructions d'utilisation	10	NL Gebruiksaanwijzing	34
ES Instrucciones de uso	14	Symbols	39
IT Istruzioni per l'uso	18	Our Teams	40
PL Sposób użycia	22		

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Influenza-Stämme

Influenzavirus Typ A	Influenzavirus Typ B
H1N1 A/NWS/33	B/Russia/69
H1N1 A/WS/33	B/Hong Kong/5/72
H1N1 A/New Jersey/8/76	B/Lee/40
H1N1 A/Mal/302/54	B/Bright
H3N2 A/Aichi/2/68	B/B/R5
H3N2 A/Hong Kong/8/68	B/Wisconsin/1/2010 (Yamagata lineage)
H3N2 A/Port Chalmers/1/73	B/Florida/78/2015 (Victoria lineage)

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14. Referenzen

- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
- https://www.rki.de/DE/Content/Infekt/EpidBull/Merkblaetter/Ratgeber_Influenza_saisonal.html
- <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4848591/>

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1. Intended Use

The NADAL® Influenza A/B Test is a lateral flow chromatographic immunoassay for the qualitative detection of influenza virus type A and B antigens (nucleoproteins) in the nasal swabs, pharyngeal swabs or nasal/nasopharyngeal aspirates of symptomatic patients. This test is intended for use as an aid in the rapid, differential diagnosis of influenza virus A and B infections and is designed for professional use only.

2. Introduction and Clinical Significance

Influenza is a highly contagious viral infection of the upper respiratory tract which is characterised by antigen variability, seasonality and the impact on the general population. Of the two main types (A and B) of influenza viruses, influenza A subtypes are differentiated by the antigen variability of surface glycoproteins (haemagglutinin and neuraminidase). The influenza A virus is the most prevalent and is associated with the most serious epidemics. Influenza can cause severe complications such as bronchitis or pneumonia, particularly in children, the elderly or in those with chronic respiratory disease. However, it most commonly occurs as a mild viral infection transmitted by respiratory secretions through sneezing or coughing. There are many other viral infections that can mimic influenza symptoms, making laboratory tests necessary to distinguish it from other acute respiratory infections.

With a sensitivity and specificity of almost 100%, RT-PCR is the gold standard of laboratory diagnostics. However, as a laboratory method, it does not allow immediate diagnosis and requires qualified personnel and expensive laboratory equipment. Antiviral influenza drugs, which have been available since the late 1990s, are most effective if administered early after the onset of the illness (within 48 hours). Rapid tests for the detection of influenza antigens can therefore enable that a successful antiviral therapy is carried out in good time, thereby improving patient healthcare and reducing health costs.

3. Test Principle

The NADAL® Influenza A/B Test is a lateral flow chromatographic immunoassay for the qualitative detection of influenza virus type A and B antigens (nucleoproteins) in nasal swabs, pharyngeal swabs or nasal/nasopharyngeal aspirates.

Antibodies specific to influenza virus type A and B nucleoproteins are pre-coated onto the test line region (A) and (B) of the membrane. During testing, the extracted specimen reacts with the antibodies to influenza virus type A and B antigens which are conjugated to coloured particles. The mixture then migrates along the membrane chromatographically by capillary action and reacts with the antibodies to influenza virus type A and B antigens in the test line region (A) and (B) of the membrane. The presence of a coloured line in the test line region (A) and/or (B) indicates a positive result. The absence of a coloured line in the test line region (A) and/or (B) indicates a negative result.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® Influenza A/B test cassettes
- Additional material provided according to 93/42/EEC:
Due to possible supply shortages of accessory medical products as a result of the current COVID-19 pandemic, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

a) 10 sterile swabs, CE0197

Miraclean Technology Co., Ltd.

No. 18, Rongshuxia Industrial Zone

Tongle Community, Longgang District Shenzhen

518116 Guangdong, China (authorised EU

representative Share Info Consultant Service LLC

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Düsseldorf, Germany)

b) 10 sterile swabs, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd

Touqiao Town, Guangling District, Yangzhou,

Jiangsu 225109 China (authorised EU

representative Llns Service & Consulting GmbH,

Obere Seegasse 34/2, 69124 Heidelberg,

Germany)

- 10 extraction tubes incl. dropper caps
- 2 control swabs: 1 Influenza A+/B- Control Swab (non-viable Influenza A, 0.02% NaN₃); 1 Influenza A-/B+ Control Swab (non-viable Influenza B, 0.02% NaN₃)
- 1 reagent holder
- 1 buffer (8 mL)
- 1 package insert

5. Additional Materials Required

- Timer
- Aspiration device with catheter and mucus trap (for aspirates)

6. Storage & Stability

Test kits should be stored at 2-30°C until the expiry date. Test cassettes are stable until the expiry date printed on the sealed foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze tests. Do not use tests beyond the expiry date indicated on the packaging.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not use specimens contaminated by blood.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.

- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

8. Specimen Collection and Preparation

Nasal swab specimen:

- Gently insert a sterile swab into the nostril.
- Collect secretion by rotating the swab several times against the turbinates.

Pharyngeal swab specimen:

- Gently insert a sterile swab into the pharynx and collect secretion by brushing the swab several times against the reddened posterior pharyngeal wall and palatine tonsils. Be careful not to collect saliva with the swab.


Nasal/nasopharyngeal aspirate specimen:

- To collect a nasal/nasopharyngeal aspirate specimen, connect a catheter to the mucus trap attached to the aspiration device. Insert the catheter through the nostril and into the nasal cavity before starting the aspiration device. Once the sample has been collected, dip a sterile swab into the nasal/nasopharyngeal aspirate. Allow the swab to absorb the specimen.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.


1. Place a clean extraction tube into the designated area of the reagent holder. Holding the buffer bottle vertically and without touching the edge of the tube, add 8 drops (approximately 300 µL) to the extraction tube.
 
2. Insert the swab with the collected specimen into the tube. Swirl the swab for approximately 10 seconds, pressing it against the wall of the extraction tube to extract the antigens contained in the swab.
 
3. Remove the swab, squeezing it firmly against the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.

4. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
5. Place the test cassette on a clean and level surface.
6. Attach a dropper cap to the tube and add 3 drops (approximately 80 µL) of the specimen to the specimen well (S) of the test cassette.
 
7. Start the timer.
8. Wait for the coloured line(s) to appear. Read the test result after 8 minutes. Do not interpret the result after more than 15 minutes.
 

10. Result Interpretation


Positive for Influenza A

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'A'. The result indicates that influenza virus type A antigens have been detected in the specimen.


Positive for Influenza B

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'B'. The result indicates that influenza virus type B antigens have been detected in the specimen.

Positive for Influenza A and Influenza B


In addition to the control line 'C', a coloured line develops in the test line region 'A' and another in the test line region 'B'. The result indicates that influenza virus type A and B antigens have been detected in the specimen.

Note: The colour intensity in the test line region (A) and (B) may vary depending on the influenza virus type A or B antigen concentration in the specimen. Therefore, any shade of colour in the test line region (A) or (B) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.


Negative

A coloured line develops in the control line region 'C'. No lines develop in the test line region 'A' and 'B'.

Invalid

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

Influenza A+/B- and Influenza A-/B+ control swabs are provided and should be tested with each test kit.

Procedure for external quality control testing:

Influenza A+/B- and Influenza A-/B+ control swabs should be tested by following the instructions described in section 9 'Test Procedure'.

If the controls do not yield the expected results, do not use the test kit and contact your distributor.

12. Limitations

- The NADAL® Influenza A/B Test is for professional *in-vitro* diagnostic use only and should only be used for the qualitative detection of influenza virus type A and B antigens in nasal swabs, pharyngeal swabs or nasal/nasopharyngeal aspirate specimens. Neither the quantitative value nor the rate of increase in influenza virus type A and B concentration can be determined with this qualitative test.
- The NADAL® Influenza A/B Test only detects the presence of viable and non-viable influenza virus type A and/or B strains in specimens.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- A negative test result obtained using the NADAL® Influenza A/B Test should be confirmed by culture. A negative test result may be obtained if the concentration of influenza virus type A and/or B in nasal swabs, pharyngeal swabs or nasal/nasopharyngeal aspirates is below the detection limit of the test.
- Excess blood or mucus in nasal swabs, pharyngeal swabs or nasal/nasopharyngeal aspirate specimens may interfere with the test performance and yield a false positive result.

- The accuracy of the test depends on the quality of the specimen. False negative results may occur due to improper specimen collection or storage.
- The use of over-the-counter and prescription nasal sprays at high concentrations can interfere with test results, leading to either invalid or incorrect test results.
- A positive result for influenza virus type A and/or B does not preclude an underlying co-infection with another pathogen; therefore the possibility of an underlying bacterial infection should be considered.
- The performance of the NADAL® Influenza A/B Test has not been evaluated for monitoring antiviral treatment for influenza.

13. Performance Characteristics

Sensitivity and specificity

The NADAL® Influenza A/B Test was evaluated using specimens obtained from patients. RT-PCR was used as a reference method. Specimens were considered positive or negative if RT-PCR showed a corresponding result.

		Influenza virus type A			influenza virus type B		
		RT-PCR		Total	RT-PCR		Total
		Positive	Negative		Positive	Negative	
Influenza A/B Test	Positive	68	14	82	49	7	56
	Negative	10	242	252	4	274	278
Total		78	256	334	53	281	334
Relative sensitivity		87.2%			92.5%		
Relative specificity		94.5%			97.5%		
Overall agreement		92.8%			96.7%		

Interfering substances

100 µL of each one-percent solution of the substances listed below were added to several swabs. Then 10 µL of a negative control, 10 µL of a strong positive or of a weak positive influenza A specimen, as well as 10 µL of a strong positive or of a weak positive influenza B specimen were added separately to these swabs. The swabs were evaluated using the NADAL® Influenza A/B Test.

Mucosolvan® Ambroxol Hydrochloride Oral Solution, Nin Jiom Pei Pa Koa (cough syrup), Flixonase™ nasal spray (fluticasone propionate), Dextromethorphan Hydrobromide Oral Solution, Hyland's 4 Kids Cold'n Cough Liquid Safe Natural Relief, Durham's Canker-Rid®, Listerine® mouthwash, Scope® mouthwash, haemoglobin, acetaminophen, Aspirin®, ibuprofen.

None of the substances interfered with the assay at the concentrations tested.

Cross-reactivity

Specimens spiked with the following pathogens were tested using the NADAL® Influenza A/B Test:

Viruses

Adenovirus, coxsackievirus, cytomegalovirus, parainfluenza virus type 1, 2, 3, 4a, enterovirus, mumps virus, respiratory syncytial virus, rhinovirus.

Bacteria

Bordetella pertussis, *Haemophilus parainfluenzae*, *Staphylococcus aureus*, *Streptococcus agalactiae*, *Neisseria meningitides*, *Streptococcus* sp. group A, B, C.

No cross-reactivity was observed with the aforementioned pathogens when tested using the NADAL® Influenza A/B Test.

Analytical reactivity with human influenza strains

The human influenza strains listed below were serially diluted to strong and weak positive concentrations. Swabs were spiked with 10 µL of each dilution, tested in triplicates using the NADAL® Influenza A/B Test and demonstrated a positive reaction.

Influenza strains

Influenza virus type A	Influenza virus type B
H1N1 A/NWS/33	B/Russia/69
H1N1 A/WS/33	B/Hong Kong/5/72
H1N1 A/New Jersey/8/76	B/Lee/40
H1N1 A/Mal/302/54	B/Brigit
H3N2 A/Aichi/2/68	B/B/R5
H3N2 A/Hong Kong/8/68	B/Wisconsin/1/2010 (Yamagata lineage)
H3N2 A/Port Chalmers/1/73	B/Florida/78/2015 (Victoria lineage)

Precision

Repeatability and reproducibility

Precision was established by testing 5 replicates of negative, influenza A as well as influenza B strong and weak positive specimens. Repeatability was established within the reproducibility study. Testing was performed using 3 independent NADAL® Influenza A/B test lots on 3 separate days.

The NADAL® Influenza A/B Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

14. References

1. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
2. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
3. CDC website: <http://www.cdc.gov/flu/>
4. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
5. https://www.rki.de/DE/Content/Infekt/EpidBull/Merkblaetter/Ratgeber_Influenza_saisonal.html
6. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4848591/>

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Souches de grippe

Virus Influenza de type A	Virus Influenza de type B
H1N1 A/NWS/33	B/Russia/69
H1N1 A/WS/33	B/Hong Kong/5/72
H1N1 A/New Jersey/8/76	B/Lee/40
H1N1 A/Mal/302/54	B/Brigit
H3N2 A/Aichi/2/68	B/B/R5
H3N2 A/Hong Kong/8/68	B/Wisconsin/1/2010 (Yamagata lineage)
H3N2 A/Port Chalmers/1/73	B/Florida/78/2015 (Victoria lineage)

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14. Bibliographie

- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
- https://www.rki.de/DE/Content/Infekt/EpidBull/Merkblaetter/Ratgeber_Influenza_saisonal.html
- <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4848591/>

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14. Referencias

1. Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
2. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
3. CDC website: <http://www.cdc.gov/flu/>
4. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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b) 10 tamponi sterili, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (rappresentante autorizzato EU Lins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)

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14. Lähteet

1. Murphy, B.R., and R.G. Webster, 1996. Orthomyxovirus, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
2. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
3. CDC website: <http://www.cdc.gov/flu/>
4. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
5. https://www.rki.de/DE/Content/Infekt/EpidBull/Merkblaetter/Ratgeber_Influenza_saisonal.html
6. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4848591/>

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nal von minden GmbH