# NADAL<sup>®</sup> SARS-CoV-2/Influenza A+B/RSV Ag Test (test cassette)

REF 243502N-20

<b>de</b> Gebrauchsanweisung	2	<b>cs</b> Návod k použití	31
<b>en</b> Instructions for use	7	<b>fi</b> Käyttöohje	35
<b>fr</b> Instructions d'utilisation	11	<b>sv</b> Användarinstruktioner	39
<b>es</b> Instrucciones de uso	15	<b>da</b> Brugervejledning	43
<b>it</b> Istruzioni per l'uso	19	<b>no</b> Bruksanvisning	47
<b>pl</b> Sposób użycia	23	<b>nl</b> Gebruiksaanwijzing	51
<b>pt</b> Instruções de Utilização	27	Symbols	55
		Our Teams	56

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#### 14. Referenzen

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.

5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
8. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
9. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C. Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
10. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
11. Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S. Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
12. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
13. Ahluwalia, G.J. Embree, P. McNicol, B. Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. J. Clin. Microbiol. 257 : 763-767.

Rev. 0, 2022-05-09 OM

### 1. Intended Use

The NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2, influenza A and B as well as respiratory syncytial virus (RSV) antigens in nasopharyngeal swab specimens from individuals with a suspected SARS-CoV-2/influenza/RSV infection in conjunction with clinical presentation and the results of other laboratory tests. The NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test is designed for *in-vitro* diagnostic and professional use only. It is intended for use in clinical laboratories and by healthcare professionals for point-of-care testing only. Not for home-use.

### 2. Introduction and Clinical Significance

COVID-19 (coronavirus disease) is a disease caused by the SARS-CoV-2 virus. It is highly contagious and spreads quickly. COVID-19 most often causes respiratory symptoms that can feel much like a cold, the flu, or pneumonia. COVID-19 may affect more than the lungs and respiratory system.

Influenza is an acute viral infection of the upper or lower respiratory tract marked by fever, chills, and a generalised feeling of weakness and pain in the muscles, along with varying degrees of soreness in the head and abdomen.

Respiratory syncytial virus, or RSV, is a common respiratory virus that usually causes mild, cold-like symptoms. Most people recover in a 1-2 weeks, but RSV can be serious. Infants and older adults are more likely to develop severe RSV symptoms and more likely to require hospitalisation. Vaccines are available to protect older adults from a severe course of RSV. Preventative options are also available to protect infants and young children from severe RSV symptoms.

### 3. Test Principle

The NADAL® SARS-CoV-2 Ag Test has one line of anti-SARS-CoV-2 antibodies in the test line region 'T' and one line of anti-mouse IgG antibodies in the control line region 'C'. When an extracted specimen is added to the specimen well, it reacts with the labeled antibodies to form complexes. The mixture then migrates along the membrane by capillary action and interacts with the precoated anti-SARS-CoV-2 antibodies in the test line region 'T'. If the specimen contains SARS-CoV-2 antigens, a red test line 'T' will appear, indicating a positive SARS-CoV-2 antigen result. Otherwise, the test result is to be considered negative. The test cassette also contains a control line 'C', which should turn red for all valid tests. If the quality control line 'C' does not appear, the test result is to be considered invalid, even if the test line 'T' appears.

The NADAL® Influenza A+B Ag Test has one line of anti-influenza A antibodies in the test line region 'A', one line of anti-influenza B antibodies in the test line region 'B' and one line of anti-mouse IgG antibodies in the control line region 'C'. When an extracted specimen is added to the specimen well, it reacts with the labeled antibodies to form complexes. The mixture then migrates along the membrane by capillary action and interacts with the precoated anti-influenza A antibodies and anti-influenza B antibodies in the corresponding test line regions. If the specimen contains influenza A or influenza B antigens, a red test line will appear, indicating the presence of influenza A or influenza B antigens. Otherwise, the test result is to be considered negative. The test cassette also contains a control

line 'C' which should turn red for all valid tests. If the control line 'C' does not appear, the test result is to be considered invalid, even if the test line appears.

The NADAL® RSV Ag Test has one line of anti-RSV antibodies in the test line region 'T' and one line of anti-mouse IgG antibodies in the control line region 'C'. When an extracted specimen is added to the specimen well, it reacts with the labeled antibodies to form complexes, the mixture then migrates along the membrane by capillary action and interacts with the coated anti-RSV antibodies in the test line region 'T'. If the specimen contains RSV antigens, a red test line will appear, indicating the presence of RSV antigens. Otherwise, the test result is to be considered negative. The test cassette also contains a control line 'C' which should turn red for all valid tests. If the control line 'C' does not appear, the test result is to be considered invalid, even if the test line appears.

### 4. Reagents and Materials Supplied

- 20 NADAL® SARS-CoV-2/Influenza A+B/RSV Ag test cassettes
- Additional material provided in accordance with 93/42/EEC: Due to possible supply shortages of accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

a) 20 sterile swabs, CE 0197

CITOTEST LABWARE MANUFACTURING CO., LTD  
No. 339 Beihai West Road, Haimen, 226100  
Jiangsu, P.R. China (authorised EU representative:  
WellKang Ltd, Enterprise Hub, NW Business  
Complex, 1 Beraghmore Rd., Derry, BT48 8SE,  
Northern Ireland)

b) 20 sterile swabs, CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room  
402, Building 2, No. 2628, Yuhangtang Road,  
Cangqian Street, Yuhang District, Hangzhou,  
311121 Zhejiang, China (authorised EU  
representative: Zoustech S.L., Paseo de la  
Castellana 141, 28049 Madrid, Spain)

- 20 extraction tubes containing extraction solution (approximately 600 µL each)\*
- 20 dropper caps
- 1 reagent holder
- 1 package insert

\*contains the following preservative: ProClin™ 300: <0.03%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) N° 1272/2008 CLP. Concentrations are below the exemption threshold of <0.03%.

### 5. Additional Materials Required

- Timer

### 6. Storage & Stability

Store the test kit in a cool, dry place between 2-30°C. Keep away from light. Exposure to temperatures and/or humidity outside the specified conditions may cause inaccurate results. Do not freeze. Use the test kit at temperatures between 15-30°C. Use the test kit between 10-90% humidity. Do not use the test kit beyond the expiration date (printed on the foil pouch and box).

**Note:** All expiration dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

## 7. Warnings and Precautions

- For *in-vitro* diagnostic use only.
- This package insert must be read completely before performing the test. Failure to follow the instructions provided in the package insert may lead to inaccurate test results.
- Do not open the sealed pouch unless ready to conduct the assay.
- Do not use expired tests.
- Bring all reagents to room temperature (15–30°C) before use.
- Do not use the components from any other type of test kits as a substitute for the components in this test kit.
- Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Dispose of all specimens and materials used to perform the test as biohazardous waste.

## 8. Specimen Collection and Preparation

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.


### Specimen Collection:

Use the nasopharyngeal swab supplied in the test kit. It is important to obtain as much secretion as possible. Therefore, to collect a nasopharyngeal swab sample, carefully insert the sterile swab into the nostril that presents the most secretions under visual inspection. Keep the swab near the floor of the nasal septum while gently pushing it into the posterior nasopharynx. Rotate the swab several times, and then remove it from the nasopharynx.


## 9. Test Procedure

**Bring tests and extraction solutions to room temperature (15–30°C) prior to testing.**

1. Open the extraction solution (in the sealed tube).
2. Collect a specimen as outlined in the section 'Specimen Collection and Preparation'.
3. Insert the swab with the collected specimen into the extraction tube filled with extraction solution. Roll the swab 5 times while pressing the tip against the bottom and wall of the extraction tube.
4. Remove the swab while squeezing the wall of the tube to extract the liquid from the swab. Try to expunge as much liquid as possible. Dispose of the used swab as biohazardous waste.
5. Attach a dropper cap to the extraction tube.

6. Take out a test cassette from its sealed foil pouch and put it on a clean and level surface.
7. Apply 3 drops of the extracted specimen into each specimen well. Please avoid bubbles forming during application.
8. Read the test result after 15 minutes. Do not read results after more than 20 minutes.

## 10. Result Interpretation

### For SARS-CoV-2:

#### Positive result:

Both the control line 'C' and the test line 'T' appear.

#### Negative result:

Only the control line 'C' appears, no line appears in the test line region 'T'.

#### Invalid result:

The control line 'C' fails to appear, indicating the test is invalid no matter whether the test line 'T' appears or not. Collect a new specimen and perform another test with a new test cassette.

### For Influenza A+B:

#### Positive results:

##### Positive for influenza A antigens:

Both the control line 'C' and the influenza A test line 'A' appear, while the influenza B test line 'B' does not appear.

##### Positive for influenza B antigens:

Both the control line 'C' and the influenza B test line 'B' appear, while the influenza A test line 'A' does not appear.


##### Positive for influenza A and B antigens:

All 3 lines appear, including the control line 'C' as well as the influenza A and influenza B test lines 'A' and 'B'.


#### Negative result:

Only the control line 'C' appears, no lines appear in the influenza A and influenza B test line regions 'A' and 'B'. This indicates the test result is negative for both influenza A and influenza B antigens.

#### Invalid results:

The control line 'C' fails to appear, indicating that the test is invalid, no matter whether the influenza A or influenza B test line appears or not.

Collect a new specimen and perform another test using a new test cassette.

#### For RSV:

##### Positive result:

Both the control line 'C' and the test line 'T' appear.

##### Negative result:

Only the control line 'C' appears, no line appears in the test line region 'T'.

##### Invalid result:

The control line 'C' fails to appear, indicating that the test is invalid, no matter whether the test line 'T' appears or not. Collect a new specimen and perform another test with a new test cassette.

### 11. Quality Control

Internal procedural controls are included in the test. A coloured line appearing in the control region 'C' is an internal procedural control. This procedural control line indicates that sufficient flow has occurred, and the functional integrity of the test cassette has been maintained. Control standards are not supplied with this test kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

### 12. Limitations

- The NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test is for *in-vitro* diagnostic use only and should not be re-used.
- The used test should be treated as potentially infectious material and be disposed of properly.
- The test kit should be kept away from direct sunlight, moisture and heat.
- Please check whether the test kit has any damage and check the expiry date before use.
- The sample volume may affect the accuracy of the test result. Inaccurate sample volume may cause a false positive or negative result.
- Test results must be evaluated in conjunction with other clinical data available to the physician.
- Positive test results do not rule out co-infections with other pathogens.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.
- Please be very careful when collecting nasopharyngeal swab specimens from children.
- Components from different batches are not to be used in combination.

### 13. Performance Characteristics

#### Clinical performance

##### Diagnostic sensitivity and specificity

A total of 362 specimens were tested using the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test. These specimens were obtained by collecting nasopharyngeal swabs from

symptomatic patients. The performance of the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test was compared to a commercialised molecular assay.

Table 1: SARS-CoV-2 performance of the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test as compared to a PCR Test.

NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test (for SARS-CoV-2)	PCR		
	Positive	Negative	Total
	Positive	102	1
Negative	4	255	259
Total	106	256	362

Diagnostic sensitivity: 96.23% (102/106, 95%CI, 90.70% - 98.52%)

Diagnostic specificity: 99.61% (255/256, 95%CI, 97.82% - 99.93%)

Overall agreement: 98.62% (357/362, 95%CI, 96.81% - 99.41%)

Table 2: Influenza A performance of the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test as compared to a PCR Test.

NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test (for influenza A)	PCR		
	Positive	Negative	Total
	Positive	55	4
Negative	3	300	303
Total	58	304	362

Diagnostic sensitivity: 94.83% (55/58, 95%CI, 85.86% - 98.23%)

Diagnostic specificity: 98.68% (300/304, 95%CI, 96.67% - 99.49%)

Overall agreement: 98.07% (355/362, 95%CI, 96.06% - 99.06%)

Table 3: Influenza B performance of the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test as compared to a PCR Test.

NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test (for Influenza B)	PCR		
	Positive	Negative	Total
	Positive	32	5
Negative	1	324	325
Total	33	329	362

Diagnostic sensitivity: 96.97% (32/33, 95%CI, 84.68% - 99.46%)

Diagnostic specificity: 98.48% (324/329, 95%CI, 96.49% - 99.35%)

Overall agreement: 98.34% (356/362, 95%CI, 96.43% - 99.24%)

Table 4: RSV performance of the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test as compared to a PCR Test.

NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test (for RSV)	PCR		
	Positive	Negative	Total
	Positive	111	2
Negative	2	247	249
Total	113	249	362

Diagnostic sensitivity:	98.23% (111/113, 95%CI, 93.78% - 99.51%)
Diagnostic specificity:	99.20% (247/249, 95%CI, 97.12% - 99.78%)
Overall agreement:	98.90% (358/362, 95%CI, 97.19% - 99.57%)

### Analytical performance

#### Detection limit

Virus type/ subtype	Concentration
SARS-CoV-2	75.5 TCID <sub>50</sub> /mL
Influenza A (H1N1)	1320 TCID <sub>50</sub> /mL
Influenza A (H3N2)	568 TCID <sub>50</sub> /mL
Influenza B (Victoria lineage)	1540 TCID <sub>50</sub> /mL
Influenza B (Yamagata lineage)	2180 TCID <sub>50</sub> /mL
RSV	93.6 TCID <sub>50</sub> /mL

#### Measuring range

No prozone effect was observed when testing SARS-CoV-2 (up to  $1.51 \times 10^5$  TCID<sub>50</sub>/mL), influenza A (up to  $5.68 \times 10^7$  TCID<sub>50</sub>/mL), influenza B (up to  $4.36 \times 10^7$  TCID<sub>50</sub>/mL) and RSV (up to  $1.17 \times 10^5$  TCID<sub>50</sub>/mL).

#### Cross-reactivity

No cross-reactivity was observed with the following, potentially cross-reactive pathogens when tested using the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test: Influenza A, influenza B, adenovirus, respiratory syncytial virus, coronavirus, MERS-Coronavirus, parainfluenza virus, rhinovirus A16, *Legionella pneumophila*, *Mycobacterium tuberculosis*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Mycoplasma pneumoniae*, *Chlamydomphila pneumoniae*, *Staphylococcus aureus*, human metapneumovirus, enterovirus, *Haemophilus influenzae*, *Candida albicans*, *Bordetella pertussis*, *Staphylococcus epidermidis*, *Pneumocystis jirovecii*, pooled human nasal wash.

The NADAL® SARS-CoV-2 Ag Test may cross-react with SARS-CoV-1 which caused the 2002-2004 SARS outbreak.

#### Interfering substances

The following substances showed no interference with the NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test:

Antiviral drugs: Zanamivir, Oseltamivir, Artemether/Lumefantrine, Doxycycline hyclate, Quinine, Lamivudine, Ribavirin, Daclatasvir.

Respiratory specimens: Mucin from bovine submaxillary glands, type I-S; EDTA-anticoagulated human blood; Biotin.

Nasal sprays or drops: Neosynephrin® (phenylephrine), Afrin® Nasal Spray (oxymetazoline), Saline nasal spray.

Homeopathic allergy relief medicine: homeopathic Zicam® Allergy Relief Nasal Gel, sodium cromoglycate, olopatadine hydrochloride.

Anti-inflammatory medication: Paracetamol (Acetaminophen), acetylsalicylic acid, ibuprofen.

Antibiotics: Mupirocin, tobramycin, erythromycin, ciprofloxacin.

#### Precision

##### Repeatability

Repeatability was established by testing 20 replicates of specimens using 3 lots of the NADAL® SARS-CoV-2/Influenza

A+B/RSV Ag Test. The NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test demonstrated acceptable repeatability.

#### Reproducibility

Reproducibility was established by testing 5 replicates of specimens. Testing was performed by 5 operators using 3 NADAL® SARS-CoV-2/Influenza A+B/RSV Ag test lots at 3 different sites on 5 days. The NADAL® SARS-CoV-2/Influenza A+B/RSV Ag Test demonstrated acceptable reproducibility.

#### 14. References

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McCung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Bibliographie

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: *Fields, Virology*, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee corza agent (CCA). II. Epidemiological aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1950 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V., Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurston, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in children. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Referencias

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leiboldt JL. Coronavirus pathogenesis. *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B, J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.W, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurston, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Bibliografia

- Cui J, Li F, Shi Lu, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016;24:490-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66: 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. J.Clin. Microbiol. 257 : 763-767.

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#### 14. Bibliografia

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol* 2016;24:90-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition. B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66: 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C. Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176: 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89: 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S. Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N. Engl. J. Med.* 308: 1443-1447.
- Taber, L.H.V. Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurston, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72: 613-618.
- Ahluwalia, G.J. Embree, P. Mc Nicol, B. Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J. Clin. Microbiol.* 257: 763-767.

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#### 14. Referências

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: *Fields, Virology*, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiological aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C. Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, E.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurston, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Reference

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: *Fields, Virology*, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C. Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B., T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V., Knight, B.E., Gilbert, H.W., McClung, S.Z., Wilson, H.J., Norton, J.M., Thurston, W.H., Gordon, R.L., Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J., Embree, P., McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Lähteet

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiology and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.
- Chanock, R.M., and L. Lindberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66: 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleve, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176: 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89: 1443-1447.
- Hall, C.B., J. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308: 1443-1447.
- Taber, L.H.V., Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72: 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. J.Clin. Microbiol. 257: 763-767.

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#### 14. Referenser

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. *Epidemiology, genetic recombination, and pathogenesis of coronaviruses*, *Trends Microbiol* 2016; 24:490-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition. B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg, 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B, J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Referencer

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/fuu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B, J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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#### 14. Referanser

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, *Epidemiology, genetic recombination, and pathogenesis of coronaviruses*, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, *Coronavirus pathogenesis*, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, *Orthomyxoviruses*, pp.1397-1445. In: *Fields, Virology*, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurston, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. Mc Nicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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

#### 14. Referenties

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol* 2016; 24:490-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res* 2011; 81:85-164.
- Murphy B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition. B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.
- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hidreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V., Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

Rev. 0, 2022-05-09 Tr. ag.

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> -diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> -diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelunumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

## Our Teams

**Germany:****Regensburg**

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

**Moers**

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

**Austria:**

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

**UK & Ireland:**

Tel: +49 941 290 10-18

Free Tel –UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

**France:**

France Tel: 0800 915 240

France Fax: 0800 909 493

**Switzerland**

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

**Belgium**

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

**Luxembourg**

Lux, Tel: 800 211 16

Lux, Fax: 800 261 79

**Spain:**

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

**Italy:**

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

**Poland:**

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

**Portugal:**

Tel: +49 941 290 10-735

Tel, Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

**Netherlands:**

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

**Nordic countries:****Denmark**

Tel: +31 703075 605

Free Tel: 808 887 53

**Finland**

Tel: +31 703075 606

Free Tel: 0800 918 263

Free Fax: 0800 918 262

**Norway**

Tel: +31 703075 605

Free Tel: 800 16 731

**Sweden**

Tel: +31 703075 605

Free Tel: 020 79 09 06

**VivaChek Biotech (Hangzhou) Co., Ltd.**

Level 2, Block 2, 146 East Chaofeng Rd.,

Yuhang Economy Development Zone,

Hangzhou, 311100, China

Email: info@vivachek.com

www.vivachek.com

Lotus NL B.V.

Koningin Julianaplein 10, 1e Verd,

2595AA, The Hague, Netherlands.

Tel: +31644168999

Email: peter@lotusnl.com

**Distributor:**

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH