NADAL[®] H. pylori Ag Test (test cassette)

REF 262902N-20

de	Gebrauchsanweisung	2	cs	Návod k použití	30
en	Instructions for use	6	fi	Käyttöohje	34
fr	Instructions d'utilisation	10	sv	Användarinstruktioner	38
es	Instrucciones de uso	14	no	Bruksanvisning	42
it	Istruzioni per l'uso	18		Symbols	47
pl	Sposób użycia	22		Our Teams	48
pt	Instruções de Utilização	26			

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5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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15. Referenzen

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.

1. Intended Use

The NADAL® H. pylori Ag Test is a lateral flow chromatographic immunoassay for the qualitative detection of *Helicobacter pylori* (*H. pylori*) antigens in human faecal specimens. The test is intended for use as an aid in the diagnosis of *H. pylori* infections (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® H. pylori Ag Test is designed for professional use only.

2. Introduction and Clinical Significance

Helicobacter pylori, also known as *Campylobacter pylori*, is a spiral-shaped, gram-negative bacterium with typical flagella. The bacterium infects the gastric mucosa by secreting the vacuolating cytotoxin A (VacA) and causes several gastro-enteric diseases such as non-ulcerous dyspepsia, active gastritis or gastric and duodenal ulcer. *H. pylori* is classified as carcinogen agent type I and can increase the risk of stomach adenocarcinoma.

Many *H. pylori* strains have been isolated. The one expressing the cytotoxin associated gene A (CagA) antigen is found to be highly immunogenic and is of utmost clinical importance because of its association with cytotoxic factors. It is widely reported in the specialist literature that infected patients possessing antibodies against CagA are at five times higher risk of gastric cancer than the reference group infected with a CagA negative bacterial strain. The presence of the CagA antibodies determines the persistence of *H. pylori* infection and ulceration.

CagA antigen associated to others, such as CagII and CagC seems to act as a starting agent of a sudden inflammatory response, which may provoke peptic ulceration, allergic episodes and decrease the efficacy of therapy.

At present several invasive and non-invasive methods are available to detect *H. pylori* infection. Invasive methods which require endoscopy of the gastric mucosa with a histologic, cultural and urease investigation are cost- and time-consuming.

Non-invasive methods include classical ELISA, immunoblotting assays and breath tests. But the latter are extremely complicated and highly unselective.

3. Test Principle

The NADAL® H. pylori Ag Test enables the detection of *H. pylori* antigens through the visual interpretation of colour development on the internal test strip. Anti-*H. pylori* antibodies are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with the anti-*H. pylori* antibodies which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is a sufficient amount of *H. pylori* antigens in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper

volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® H. pylori Ag test cassettes
- 20 specimen collection tubes containing buffer*
- 1 package insert

*containing the following preservative: ProClin™ 300: <0.03%

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below the exemption threshold.

5. Additional Materials Required

- Specimen collection containers (suitable for specimen material to be tested)
- Optional: stool collectors
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.

- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed in accordance with local regulations.

8. Specimen Collection and Preparation

Notes on specimen collection and storage:


The NADAL® H. pylori Ag Test is intended for use with human faecal specimens only. Faecal specimens should be collected in clean and dry specimen collection containers containing no preservatives or transport media. Urine and the excessive dilution of specimens with toilet water can lead to inaccurate test results.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Diluted and undiluted faecal specimens can be stored at 15-30°C for up to 1 hour and at 2-8°C for up to 72 hours. Bring specimens to room temperature prior to testing.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Specimen preparation:

1. Hold the specimen collection tube upright and remove the applicator stick by unscrewing the blue cap. Be careful not to spill or splash buffer solution from the collection tube.
2. Collect specimens by inserting the applicator stick into 3 different areas of the faecal specimen to collect. Avoid scooping the faecal specimen.
Note: Please be sure to insert the applicator stick 3 times in a row into the faecal specimen. Do NOT put the applicator stick back into the tube in between collecting the 3 specimens. Be careful not to spill fluid from the tube. The quality of the specimen and compliance with these instructions may affect your test result.
3. Put the applicator stick, along with the faecal specimen, back into the specimen collection tube and tightly screw the blue cap closed.
4. Shake the specimen collection tube to completely mix the specimen and buffer.

1. Remove the NADAL® H. pylori Ag test cassette from the sealed foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

2. Place the test cassette on a clean and level surface.


3. Shake the specimen collection tube to ensure good specimen dispersion.

4. Unscrew the transparent cap. Holding the specimen collection tube vertically, transfer 3 drops of the specimen solution into the specimen well (S) of the test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

5. Start the timer.


6. Wait for the coloured line(s) to appear. Read the test result after 10 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

9. Test Procedure

Bring tests, specimens extracted in buffer and/or controls to room temperature (15-30°C) prior to testing.

11. Quality Control

The internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® H. pylori Ag Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of *H. pylori* antigens in human faecal specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of *H. pylori* antigens can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- The NADAL® H. pylori Ag Test only detects the presence of *H. pylori* antigens in specimens and should not be used as the sole criterion for a diagnosis.
- During certain antibiotic treatments, the concentration of *H. pylori* antigens may decrease below the detection limit of the test. Therefore, a diagnosis should be made with caution during antibiotic treatment.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of an *H. pylori* infection.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® H. pylori Ag Test was evaluated using clinical faecal specimens in comparison to such reference methods as biopsy (bacterial culture), histology and a rapid urease test (RUT).

The results are presented in the following table:

NADAL® H. pylori Ag Test	Biopsy/Histology/RUT			
		Positive	Negative	Total
	Positive	114	2	116
Negative	2	195	197	
Total	116	197	313	

Diagnostic sensitivity: 98.28% (93.91% - 99.79%)*

Diagnostic specificity: 98.98% (96.38% - 99.88%)*

Overall agreement: 98.72% (96.76% - 99.65%)*

*95% confidence interval

Analytical performance

Detection limit

The detection limit of the NADAL® H. pylori Ag Test is 1 ng/mL recombinant *H. pylori* antigen.

Measuring range

No adverse effect on T-line formation (prozone effect) was observed when testing specimens containing an antigen concentration as high as 500 ng/mL recombinant *H. pylori* antigen. Thus, the measuring range of the test is between 1 ng/mL and at least 500 ng/mL. Higher concentrations have not been tested. If you suspect a prozone effect, dilute the specimen with the buffer from a new collection tube and repeat the measurement.

Analytical specificity

Cross-reactivity study

Negative specimens were spiked with the following, potentially cross-reacting bacteria and viruses and tested at the concentration of 1.0×10^9 microorganisms/mL using the NADAL® H. pylori Ag Test:

<i>Staphylococcus aureus</i>	<i>Pseudomonas aeruginosa</i>	<i>Enterococcus faecalis</i>
Group C <i>Streptococcus</i>	<i>Klebsiella pneumoniae</i>	<i>Moraxella catarrhalis</i>
<i>Candida albicans</i>	<i>Proteus mirabilis</i>	<i>Acinetobacter spp</i>
<i>Salmonella enterica</i>	<i>Gardnerella vaginalis</i>	<i>Acinetobacter calcoaceticus</i>
<i>E. coli</i>	<i>Chlamydia trachomatis</i>	<i>Neisseria gonorrhoeae</i>
Group B <i>Streptococcus</i>	<i>Proteus vulgaris</i>	<i>Enterococcus faecium</i>
<i>Haemophilus influenzae</i>	<i>Neisseria meningitidis</i>	<i>Rotavirus</i>

No cross-reactivity with the specimens was observed when tested using the NADAL® H. pylori Ag Test.

Interference study

Negative and low positive specimens spiked with the following potentially interfering substances at the concentrations listed below showed no interference with the NADAL® H. pylori Ag Test.

Substance	Concentration	Substance	Concentration
Haemoglobin	10 mg/mL	Ibuprofen	250 µg/mL
Acetaminophen	0.2 mg/mL	Naproxen	250 µg/mL
Acetylsalicylic acid	0.2 mg/mL	Salicylic acid	250 µg/mL
Bilirubin	10 mg/mL	Metronidazole	0.12 mg/mL
Triglycerides	50 mg/mL	Chloroquine	0.2 mg/mL
Human serum albumin	200 mg/mL	Efavirenz	0.12 mg/mL
Creatine	1 mg/mL	Rifampicin	50 µg/mL
Cyclobenzaprine	250 µg/mL	Ethanol	1%

Precision**Repeatability**

Repeatability was established by testing 20 replicates of negative, low, medium and high positive specimens using 3 independent lots of the NADAL® H. pylori Ag tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5 - 100%). The NADAL® H. pylori Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive specimens. Testing was performed by 5 operators using 3 independent NADAL® H. pylori Ag test lots on 5 separate days. >99% of the specimens were correctly identified (375/375 correct tests per concentration, 95% confidence interval: 99% - 100%). The NADAL® H. pylori Ag Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the event of any serious incidents related to the performance of the NADAL® H. pylori Ag Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. *Med. J. Australia.* (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. *New England J. Med.* (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. *Amer. J. Gastroenterology.* (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. *Am j. Med.* 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. *Am J Gastroenterol.* 1996;91:1112-1115.

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15. Bibliographie

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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15. Referencias

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-41S.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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15. Bibliografia

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am J. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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15. Lähteet

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. *Med. J. Australia*, (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. *New England J. Med.* (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. *Amer. J. Gastroenterology*. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. *Am j. Med.* 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. *Am J Gastroenterol.* 1996;91:1112-1115.

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