nal von minden Drug-Screen[®]

Saliva Tracer

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Rev.1.03; 2025-01-22 MaKa

1. Intended Use

The nal von minden Drug-Screen® Saliva Tracer is a lateral flow chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in saliva at the following cut-off concentrations:

Test	Calibrator	Cut-off (ng/mL)
6-MAM (6-Monoacetylmorphine)	6-Monoacetylmorphine	10
ABP (AB-PINACA)	AB-PINACA pentanoic acid metabolite	10
AMP (Amphetamine)	d-Amphetamine	50
BAR (Barbiturates)	Secobarbital	50
BUP (Buprenorphine)	Buprenorphine	5/10
BZD (Benzodiazepines)	Oxazepam	10/20
COC (Cocaine)	Benzoylecgonine	10/20/50
COT (Cotinine)	Cotinine	30
FYL (Fentanyl)	Fentanyl	10
KET (Ketamine)	Ketamine	30/50
MDMA (Methylenedioxy-methamphetamine)	d,l-Methylenedioxy-methamphetamine	50
MET (Methamphetamine)	d-Methamphetamine	50
MOR/OPI (Opiates)	Morphine	10/40/50
MTD (Methadone)	Methadone	30
OXY (Oxycodone)	Oxycodone	20/40
PCP (Phencyclidine)	Phencyclidine	10
SPC (Synthetic Cannabinoids)	JWH-018 5-Pentanoic acid metabolite	25
THC (Cannabis)	Δ ⁹ -THC	15/40
TML (Tramadol)	Cis-Tramadol	30

This assay provides only a preliminary analytical test result. A more specific alternate chemical method should be used to confirm a preliminary positive analytical result. Gas chromatography/mass spectrometry (GC/MS), gas chromatography/tandem mass spectrometry (GC/MS/MS), liquid chromatography/mass spectrometry (LC/MS) or liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods. Professional judgment should be applied to any drug of abuse screen test result, particularly when preliminary positive results are indicated.

2. Summary

The nal von minden Drug-Screen® Saliva Tracer is a rapid saliva screening test that can be performed without the use of a device. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in human saliva.

6-Monoacetylmorphine (6-MAM)

6-Monoacetylmorphine (6-MAM) or 6-Acetylmorphine (6-AM) is one of three active metabolites of heroin (diacetylmorphine), the others being morphine and the much less active 3-Monoacetylmorphine (3-MAM). 6-MAM occurs as a metabolite of heroin, which is rapidly created from heroin in the body. Heroin is rapidly metabolized by esterase enzymes in the brain and has an extremely short half-life. It has also relatively weak affinity to μ -opioid receptors because the 3-hydroxy group, essential for effective binding to the receptor, is masked by the acetyl group. Therefore, heroin acts as a pro-drug, serving as a lipophilic transporter for the systemic delivery of morphine, which actively binds with μ -opioid receptors.

AB-PINACA (ABP)

Synthetic cannabinoids are designer drugs that are structurally different from THC (the active component of cannabis) but act in similar ways to affect the cannabinoid receptor system in the brain. Over the past few years, this class of designer drugs has

mainstreamed to become globally popular and increasingly problematic. Synthetic cannabinoids fall into seven major structural groups:

1. Naphthoylindoles (e.g. JWH-018, JWH-073)
2. Naphthylmethylnindoles (JWH-175, JWH-184, JWH-185, JWH-199)
3. Naphthoylpyrroles (JWH-145, JWH-146, JWH-147, etc)
4. Naphthylmethylnindenes (JWH-176)
5. Phenylacetylindoles (JWH-250, JWH-251, JWH-302)
6. Cyclohexylphenols (e.g. CP-47,497)
7. Dibenzopyrans (classic cannabinoid structure such as HU-210 and HU-211)

New structural group: Aminoalkylindazoles (AB-PINACA, AB-FUBINACA, AB-CHMINACA, etc)

In their original, chemical state, synthetic cannabinoids are liquid. The drugs are usually sold combined with dried herbs that emulate marijuana and are intended for smoking although powdered versions are also available. As laws are written to control these drugs with each new synthetic cannabinoid class as they are introduced to the market, the older versions (JWH-018, JWH-073) are seen less frequently than years past. The current trend shows the aminoalkylindazole based drugs such as AB-PINACA, AB-FUBINACA and AB-CHMINACA.

Amphetamine (AMP)

Amphetamine is a sympathomimetic amine with therapeutic indications. The drug is often self-administered by nasal inhalation or oral ingestion. Depending on the route of administration, amphetamine can be detected in oral fluid as early as 5-10 minutes following use¹. Amphetamine can be detected in oral fluids for up to 72 hours after use¹.

Barbiturates (BAR)

Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of Barbiturates leads to tolerance and physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death.

Buprenorphine (BUP)

Buprenorphine is a potent analgesic often used in the treatment of opioid addiction. The drug is sold under the trade names Subutex™, Buprenex™, Temgesic™, and Suboxone™ which contain Buprenorphine HCl alone or in combination with Naloxone HCl. Therapeutically, Buprenorphine is used as a substitution treatment for opioid addicts. Substitution treatment is a form of medical care offered to opiate addicts (primarily heroin addicts) based on a similar or identical substance to the drug normally used. In substitution therapy, Buprenorphine is as effective as Methadone but demonstrates a lower level of physical dependence.

Substantial abuse of Buprenorphine has also been reported in many countries where various forms of the drug are available. The drug has been diverted from legitimate channels through theft, doctor shopping and fraudulent prescriptions, and been abused via intravenous, sublingual, intranasal and inhalation routes.

Benzodiazepines (BZD)

Benzodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep disorders. They produce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Because they are safer and more effective, Benzodiazepines have replaced Barbiturates in the treatment of both anxiety and insomnia. Benzodiazepines are also used as sedatives before some surgical and medical procedures, and for the treatment of seizure disorders and alcohol withdrawal. Risk of physical dependence increases if Benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal discomfort, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception.

Cocaine (COC)

Cocaine is a potent central nervous system (CNS) stimulant and a local anesthetic derived from the coca plant (erythroxylum coca). The drug is often self-administered by nasal inhalation, intravenous injection and free-base smoking. Depending on the route of administration, cocaine and its metabolites benzoylecgonine and ecgonine methyl ester can be detected in oral fluid as early as 5-10 minutes following use¹. Cocaine and benzoylecgonine can be detected in oral fluids for up to 24 hours after use¹.

Cotinine (COT)

Cotinine is the first-stage metabolite of nicotine, a toxic alkaloid that produces stimulation of the autonomic ganglia and central nervous system when in humans. Nicotine is a drug to which virtually every member of a tobacco-smoking society is exposed whether through direct contact or second-hand inhalation. In addition to tobacco, nicotine is also commercially available as the active ingredient in smoking replacement therapies such as nicotine gum, transdermal patches and nasal sprays. Although nicotine is excreted in saliva, the relatively short half-life of the drug makes it an unreliable marker for tobacco use. Cotinine, however, demonstrates a substantially longer half-life than nicotine, bears a high correlation with plasma cotinine levels and has been found to be the best marker for smoking status compared with saliva nicotine measurement, breath carbon monoxide testing and plasma thiocyanate testing. The window of detection for cotinine in saliva at a cut off level of 30 ng/mL is expected to be up to 1-2 days after nicotine use.

Fentanyl (FYL)

Fentanyl belongs to powerful narcotic analgesics and is a special μ -opioid receptor stimulant. Fentanyl is one of the varieties that has been listed in management of United Nations "Single Convention of narcotic drug in 1961". Among the opiate agents under international control, fentanyl is one of the most commonly used to manage moderate to severe pain. After continuous injection of fentanyl, the affected person will suffer from protracted opioid abstinence syndrome, such as ataxia and irritability etc, which illustrates the addiction after taking fentanyl in a long time. Compared with amphetamine addicts, drug addicts who take fentanyl mainly have got the possibility of a higher infection rate of HIV, show more dangerous injection behaviour and higher risk of lifelong medication overdose.

Ketamine (KET)

Ketamine is a dissociative anesthetic developed in 1963 to replace PCP (Phencyclidine). While Ketamine is still used in human anaesthesia and veterinary medicine, it is becoming increasingly abused as a streetdrug. Ketamine is molecularly similar to PCP and thus creates similar effects including numbness, loss of coordination, a sense of invulnerability, muscle rigidity, aggressive / violent behaviour, slurred or blocked speech, exaggerated sense of strength, and a blank stare. There is depression of respiratory function but not of the central nervous system, and cardiovascular function is maintained. The effects of Ketamine generally last 4-6 hours following use.

Methylenedioxymethamphetamine (MDMA)

Methylenedioxymethamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlander, 1990). The most pervasive effect of MDMA, occurring in virtually all people who take a sufficient dose of the drug, is clenching of the jaws.

Methamphetamine (MET)

Methamphetamine is a potent stimulant chemically related to amphetamine but with greater CNS stimulation properties. The drug is often self-administered by nasal inhalation, smoking or oral ingestion. Depending on the route of administration, methamphetamine can be detected in oral fluid as early as 5-10 minutes following use¹. Methamphetamine can be detected in oral fluids for up to 72 hours after use¹.

Opiates (MOR/OPI)

The drug class opiates refers to any drug that is derived from the opium poppy, including naturally occurring compounds such as morphine and codeine and semi-synthetic drugs such as heroin. Opiates act to control pain by depressing the central nervous system. The drugs demonstrate addictive properties when used for sustained periods of time; symptoms of withdrawal may include sweating, shaking, nausea and irritability. Opiates can be taken orally or by injection routes including intravenous, intramuscular and subcutaneous; illegal users may also take them intravenously or by nasal inhalation. Using an immunoassay cut-off level, codeine can be detected in oral fluid within 1 hour following a single oral dose and can remain detectable for 7-21 hours after the dose². Morphine is found more prevalently in excreted unmetabolized and is also the major metabolic product of codeine and heroin.

Methadone (MTD)

Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiate dependence (heroin, Vicodin, Percocet, morphine). Methadone is a long-acting pain reliever producing effects that last

from twelve to forty-eight hours. Ideally, methadone frees the affected person from the pressures of obtaining illegal heroin, from the dangers of injection, and from the emotional roller coaster that most opiates produce. Methadone, if taken for long periods and at large doses, can lead to a very long withdrawal period. The withdrawal symptoms from methadone are more prolonged and troublesome than those provoked by heroin cessation, yet the substitution and phased removal of methadone is an acceptable method of detoxification for patients and therapists.

Oxycodone (OXY)

Oxycodone is a semi-synthetic opioid with a structural similarity to codeine. The drug is manufactured by modifying thebaine, an alkaloid found in the opium poppy. Oxycodone, like all opiate agonists, provides pain relief by acting on opioid receptors in the spinal cord, brain, and possibly directly in the affected tissues. Oxycodone is prescribed for the relief of moderate to high pain under the well-known pharmaceutical trade names of OxyContin®, Tylox®, Percodan® and Percocet®. While Tylox®, Percodan® and Percocet® contain only small doses of oxycodone hydrochloride combined with other analgesics such as acetaminophen or aspirin, OxyContin® consists solely of oxycodone hydrochloride in a time-release form. Oxycodone is known to metabolize by demethylation into oxymorphone and noroxycodone.

Phencyclidine (PCP)

Phencyclidine, the hallucinogen commonly referred to as Angel Dust, can be detected in saliva as a result of the exchange of the drug between the circulatory system and the oral cavity. In a paired serum and saliva sample collection of 100 patients in an Emergency Department, PCP was detected in the saliva of 79 patients at levels as low as 2 ng/mL and as high as 600 ng/mL⁴.

Synthetic Cannabinoids (SPC)

Synthetic cannabinoids or K2 is a psychoactive herbal and chemical product that, when consumed, mimics the effects of cannabis. It is best known by the brand names K2 and Spice, both of which have largely become genericized trademarks used to refer to any synthetic Marijuana product. The studies suggest that synthetic marijuana intoxication is associated with acute psychosis, worsening of previously stable psychotic disorders, and also may have the ability to trigger a chronic (long-term) psychotic disorder among vulnerable individuals such as those with a family history of mental illness.

As of March 1, 2011, five cannabinoids, JWH-018, JWH-073, CP-47, JWH-200 and cannabicyclohexanol, are now illegal in the USA because these substances have the potential to be extremely harmful and therefore pose an imminent hazard to public safety.

Cannabis (THC)

THC (Δ^9 -tetrahydrocannabinol) is the primary active ingredient in cannabis (marijuana). When smoked or orally administered, THC produces euphoric effects. Users have impaired short-term memory and impaired learning ability. They may also experience transient episodes of confusion and anxiety. Long-term, relatively heavy use may be associated with behavioural disorders. The parental THC, also known as Δ^9 -THC, is present in oral

fluid after use. The detection of the drug is thought to be primarily due to the direct exposure of the drug to the mouth (oral and smoking administrations) and the subsequent sequestering of the drug in the buccal cavity³. Historical studies have shown a window of detection for THC in saliva of up to 14 hours after drug use³.

Tramadol (TML)

Tramadol (TML) is a quasi-narcotic analgesic used in the treatment of moderate to severe pain. It is a synthetic analog of codeine but has a low binding affinity to the μ -opioid receptors. Large doses of tramadol can develop tolerance and physiological dependency and lead to its abuse. Tramadol is extensively metabolized after oral administration. The major pathways appear to be N- and O-demethylation, glucuronidation or sulfation in the liver.

3. Assay Principle

The nal von minden Drug-Screen® Saliva Tracer is an immunoassay based on the principle of competitive binding. Drugs that may be present in the oral fluid specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a portion of the oral fluid specimen migrates upward by capillary action. A drug, if present in the oral fluid specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible coloured line will form in the test line region of the specific drug strip. The presence of a drug above the cut-off concentration in the oral fluid specimen will saturate all the binding sites of the antibody. Therefore, the coloured line will not form in the test line region. A drug-positive oral fluid specimen will not generate a coloured line in the specific test line region of the strip because of drug competition, while a drug-negative oral fluid specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a coloured line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents

Each test line contains anti-drug antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polyclonal antibodies and rabbit IgG.

5. Materials

Materials Provided:

- Test cassettes
- Package insert
- Procedure card

Materials Required but Not provided:

- Timer

6. Precautions

- Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- Saliva is not classified as biological hazard unless derived from a dental procedure.
- The used collector and cassette should be disposed of according to federal, state and local regulations.

7. Storage and Stability

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test cassettes must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

8. Specimen Collection and Preparation

The oral fluid specimen should be collected using the collector provided with the kit. Follow the detailed Directions for Use below. No other specimen collector should be used with this assay. Oral fluid collected at any time of the day may be used.


9. Directions for Use

Allow the test cassette, specimen, and/or controls to reach room temperature (15-30°C) prior to testing. Instruct the test subject to not place anything in the mouth including food, drink, gum or tobacco products for at least 10 minutes prior to collection.


1. Bring the pouch to room temperature before opening. Remove the test from the sealed pouch and use within one hour of opening.
2. Instruct the test subject to place their tongue against the base of the upper or lower jaw and collect saliva in the mouth.
3. Remove the collector from the cassette, then remove the cap from the collector.

4. Instruct the test subject to place the sponge of the collector between the lower cheek and gum and gently rub back and forth between the left and right cheeks and the gum until the sponge is completely saturated with saliva. Do not bite, suck, or chew the sponge as it may break.


5. Remove the collector from the mouth when two red/pink lines appear on the back of the mouth of the collector or when the red/pink lines cover three sides of the square (approx. 30 seconds). If the saturation indicator has turned red, insert the collector into the cassette. If the saturation indicator has


not turned red, place the sponge back in the mouth and continue to collect saliva until the indicator turns red.


Note: When inserting the collector into the cassette, first insert the protruding part of the collector head into the hole reserved at the sampling site (a), and then press down the handle of the collector to secure it (b).

6. Slide the cover in the direction of the arrow until the cover is blocked.

7. Place the device on a flat surface while the test is running. Negative results can be read as soon as visible lines form in both the C and T zones of the test. Read presumptive positive results after **10 minutes. Do not read results after more than 20 minutes.**

10. Interpretation of Results


Negative*: Two lines appear. One coloured line should be in the control region (C), and another adjacent coloured line should be in the test region (Drug/T). This negative result indicates that the drug concentration is below the detectable level.

*NOTE: The shade of colour in the test line region (Drug/T) will vary, but it should be considered negative whenever there is even a faint line.

Positive: One coloured line appears in the control region (C).

No line appears in the test region (Drug/T). This positive result indicates that the drug concentration is above the detectable level.

Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact the manufacturer.

11. Quality Control

A procedural control is included in the test. A coloured line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations of the Test

1. The nal von minden Drug-Screen® Saliva Tracer provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmatory result. Gas chromatography/mass spectrometry (GC/MS), gas chromatography/tandem mass spectrometry (GC/MS/MS), liquid chromatography/mass spectrometry (LC/MS) or liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods.
2. A positive test result does not indicate the concentration of drugs in the specimen or the route of administration.
3. A negative result may not necessarily indicate a drug-free specimen. Drugs may be present in the specimen below the cut-off level of the assay.

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Performance Characteristics

Accuracy

Assemble each single test into the cassette before testing, and evaluate the cassette with approximately 44-280 specimens per drug type previously collected from subjects presenting for drug screen testing which were confirmed by GC/MS. These specimens were randomized and tested using the nal von minden Drug-Screen® Saliva Tracer. Specimens were interpreted as either positive or negative at 10 minutes. The test results are shown in table below.

Table: Specimen Correlation

Method		GC/MS		% agreement with GC/MS	% total agreement with GC/MS
Saliva Tracer Test		Positive	Negative		
6-MAM 10	Positive	36	0	>99%	>99.0%
	Negative	0	128	>99%	
ABP10	Positive	4	0	>99%	>99.0%
	Negative	0	40	>99%	
AMP50	Positive	90	6	94.7%	94.8%
	Negative	5	109	94.8%	
BAR50	Positive	80	6	96.4%	95.7%
	Negative	3	121	95.3%	
BUP10	Positive	86	5	95.6%	95.7%
	Negative	4	115	95.8%	
BUP5	Positive	86	5	95.6%	95.7%
	Negative	4	115	95.8%	
BZD20	Positive	94	5	94.0%	94.8%
	Negative	6	105	95.5%	
BZD10	Positive	94	5	94.0%	94.8%
	Negative	6	105	95.5%	
COC50	Positive	38	2	95.0%	96.7%
	Negative	3	107	97.3%	
COC20	Positive	38	2	95.0%	96.7%
	Negative	3	107	97.3%	
COC10	Positive	37	3	95.0%	95.3%
	Negative	4	106	97.3%	
COT30	Positive	131	2	99.2%	98.7%
	Negative	1	96	98.0%	
FYL10	Positive	53	1	98.1%	96.7%
	Negative	4	92	95.8%	
KET50	Positive	90	6	93.8%	94.8%
	Negative	5	109	95.6%	
KET 30	Positive	49	3	94.2%	94.5%
	Negative	5	88	94.6%	
MDMA 50	Positive	96	1	97.0%	98.3%
	Negative	3	130	99.2%	
MET50	Positive	126	4	99.2%	98.2%
	Negative	1	149	97.4%	
MOR/OPI 50	Positive	89	7	93.7%	93.8%
	Negative	6	108	93.9%	
MOR/OPI 40	Positive	89	7	93.7%	93.8%
	Negative	6	108	93.9%	
MOR/OPI 10	Positive	88	8	92.6%	92.9%
	Negative	7	107	93.0%	
MTD30	Positive	116	3	97.5%	97.4%
	Negative	3	108	97.3%	
OXY40	Positive	93	0	>99%	>99.0%
	Negative	0	137	>99%	
OXY20	Positive	91	1	97.8%	98.7%
	Negative	2	136	99.3%	
PCP10	Positive	107	2	96.4%	97.4%
	Negative	4	117	98.3%	
SPC25	Positive	52	2	96.3%	96.0%
	Negative	4	92	95.8%	
THC40	Positive	45	0	95.7%	98.0%
	Negative	2	52	>99%	
THC15	Positive	43	0	95.6%	97.8%
	Negative	2	45	99%	
TML30	Positive	89	0	>99%	>99.0%
	Negative	0	121	>99%	

Analytical Sensitivity

A Phosphate-buffered saline (PBS) pool was spiked with drugs to target concentrations of $\pm 50\%$ cut-off, $\pm 25\%$ cut-off and $+300\%$ cut-off and tested with the nal von minden Drug-Screen® Saliva Tracer. The results are summarized below.

Drug conc. (Cut-off range)	n	6-MAM10		ABP10		AMP50		BAR50	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	27	3	27	3	27	3	23	7
Cut-off	30	14	16	15	15	15	15	16	14
+25% Cut-off	30	4	24	3	27	7	23	6	24
+50% Cut-off	30	0	30	0	30	0	30	0	30
+300 % Cut-off	30	0	30	0	30	0	30	0	30

Drug conc. (Cut-off range)	n	BUP10		BUP5		BZD20		BZD10	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	27	3	27	3	26	4	27	3
Cut-off	30	15	15	15	15	14	16	15	15
+25% Cut-off	30	7	23	7	23	5	25	7	23
+50% Cut-off	30	0	30	0	30	0	30	0	30
+300 % Cut-off	30	0	30	0	30	0	30	0	30

Drug conc. (Cut-off range)	n	COC50		COC20		COC10		COT30	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	25	5	25	5	26	4	25	5
Cut-off	30	15	15	15	15	15	15	20	10
+25% Cut-off	30	3	27	3	27	3	27	7	23
+50% Cut-off	30	0	30	0	30	0	30	0	30
+300 % Cut-off	30	0	30	0	30	0	30	0	30

Drug conc. (Cut-off range)	n	FYL10		KET50		KET30		MDMA50	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	24	6	25	5	25	5	26	4
Cut-off	30	15	15	18	12	16	14	19	11
+25% Cut-off	30	3	27	8	22	4	26	6	24
+50% Cut-off	30	0	30	0	30	0	30	0	30
+300 % Cut-off	30	0	30	0	30	0	30	0	30

Drug conc. (Cut-off range)	n	MET50		MOR50		MOR40		MOR10	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	28	2	27	3	27	3	26	4
Cut-off	30	16	14	15	15	13	17	13	17
+25% Cut-off	30	6	24	8	22	7	23	7	23
+50% Cut-off	30	0	30	0	30	0	30	0	30
+300 % Cut-off	30	0	30	0	30	0	30	0	30

Drug conc. (Cut-off range)	n	MTD30		OXY40		OXY20		PCP10	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	25	5	25	3	27	3	25	5
Cut-off	30	15	15	15	15	20	10	14	16
+25% Cut-off	30	7	23	7	22	4	26	10	20
+50% Cut-off	30	16	30	0	30	0	30	0	30
+300 % Cut-off	30	6	30	0	30	0	30	0	30

Drug conc. (Cut-off range)	n	SPC25		THC40		THC15		TML30	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	26	4	26	4	27	3	25	5
Cut-off	30	15	15	12	18	12	18	14	16
+25% Cut-off	30	4	26	5	25	5	25	4	26
+50% Cut-off	30	0	30	0	30	0	30	0	30
+300 % Cut-off	30	0	30	0	30	0	30	0	30

Analytical Specificity

The following table lists the concentration of compounds (ng/mL) above which the nal von minden Drug-Screen® Saliva Tracer identified positive results at a reading time of 10 minutes.

Compound	ng/mL
6-MONOACETYLMORPHINE (6-MAM10)	
6-Monoacetylmorphine	10
Diacetylmorphine(herion)	25
AB-Pinaca (ABP10)	
AB-PINACA pentanoic acid metabolite	10
ADB-PINACA N-(4-hydroxypentyl) metabolite	15
5-fluoro AB-PINACA N-(4-hydroxypentyl)	20
AB-PINACA N-(5-hydroxypentyl) metabolite	30
AB-PINACA	100
5-fluoro ADB-PINACA	250
AB-PINACA N-(4-hydroxypentyl) metabolite	10
ADB-PINACA N-(5-hydroxypentyl) metabolite	20
ADB-PINACA pentanoic acid metabolite	20
5-fluoro AB-PINACA	50
AB-FUBINACA	150
5-chloro AB-PINACA	1000

AMPHETAMINE (AMP50)	
d-Amphetamine	50
d/l-Amphetamine	100
β-Phenylethylamine	25,000
Tryptamine	12,500
p-Hydroxyamphetamine	100
(+)-3,4-Methylenedioxyampheta mine (MDA)	100
l-Amphetamine	25,000
Methoxyphenamine	12,500
BARBITURATES (BAR50)	
Amobarbital	500
5,5-Diphenylhydantoin	1000
Allobarbital	75
Barbital	1000
Talbutal	5
Butalbital	1000
Phenobarbital	50
Cyclopentobarbital	4170
Pentobarbital	1000
Alphenol	50
Aprobarbital	75
Butabarbital	25
Butethal	75
Secobarbital	50
BUPRENORPHINE (BUP10)	
Norbuprenorphine	180
Buprenorphine-3-β-D-glucuronide	100
Buprenorphine	10
Norbuprenorphine-3-β-D-glucu ronide	600
BUPRENORPHINE (BUP5)	
Norbuprenorphine	90
Buprenorphine-3-β-D-glucuronide	50
Buprenorphine	5
Norbuprenorphine-3-β-D-glucu ronide	300
BENZODIAZEPINES (BZD20)	
Oxazepam	20
Alprazolam	200
Chlordiazepoxide	100
Desalkylflurazepam	1,000
Estazolam	160
Furosemide	10,000
Midazolam	2,000
Nefopam	2,000
Norchlordiazepoxide	50
Pheniramine	100,000
α-Hydroxyalprazolam	100
7-Amino-clonazepam	10,000
Bromazepam	20
Clonazepam	2,000
Diazepam	100
Flunitrazepam	1,000
Lorazepam	1,400
Midazolam Maleate	5,000
Nitrazepam	50
Oxolinic acid	100,000
Theophylline	100,000
BENZODIAZEPINES (BZD10)	
Oxazepam	10
Alprazolam	100
Chlordiazepoxide	50
Desalkylflurazepam	500
Estazolam	80
Furosemide	5,000
Midazolam	1,000
Nefopam	1,000
Norchlordiazepoxide	25
Pheniramine	50,000
α-Hydroxyalprazolam	50
7-Amino-clonazepam	5,000
Bromazepam	10
Clonazepam	1,000
Diazepam	50
Flunitrazepam	500
Lorazepam	700
Midazolam Maleate	2,500
Nitrazepam	25
Oxolinic acid	50,000
Theophylline	50,000

COCAINE (COC50)	
Cocaine HCl	50
Benzoylcegonine	50
EcgonineHCl	37,5
Cocaeethylene	75
COCAINE (COC20)	
Cocaine HCl	20
Benzoylcegonine	20
EcgonineHCl	15
Cocaeethylene	30
COCAINE (COC10)	
Cocaine HCl	10
Benzoylcegonine	10
EcgonineHCl	7,5
Cocaeethylene	15
COTININE (COT30)	
(-) Cotinine	30
(-) Nicotine	450
FENTANYL (FYL10)	
Fentanyl	10
Perphenazine	20,000
Norfentanyl	4
KETAMINE (KET50)	
Ketamine	50
Tetrahydrozoline	20
Benzphetamine	1250
d-Methamphetamine	1250
(+) Chlorpheniramine	1250
l-Methamphetamine	2500
Clonidine	5000
Methoxyphenamine	625
Disopyramide	625
d-Norpropoxyphene	625
EDDP	2500
Pentazocine	1250
Mephentermine	1250
Phencyclidine	625
(1R, 2S) - (-)-Ephedrine	5000
Promazine	1250
4-Hydroxyphencyclidine	2500
Promethazine	1250
Levorphanol	2500
Thioridazine	2500
MDE	2500
Meperidine	1250
Dextromethorphan	75
(+)-3,4-Methylenedioxyamphetamine (MDMA)	5000
KETAMINE (KET30)	
Ketamine (KET)	30
(±)-Chlorpheniramine	50,000
Levorphanol	50
Meperidine (Pethidine)	50,000
Naloxone	10,000
Naltrexone	2,500
EDDP [2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine]	5,000
Normorphine	50,000
Oxymorphone	1,000
Pheniramine	50,000
Norketamine	400
Pantoprazole Sodium	50,000
hydromorphone	2,500
Promethazine	50,000
d-Pseudoephedrine	100,000
Phencyclidine	100
Tetrahydrozoline	5,000
Heroin (diacetylmorphine)	50,000
Methamphetamine Hydrochloride	50,000
R (-)-Methamphetamine	50,000
METHYLENEDIOXYMETHAMPHETAMINE (MDMA50)	
(±) 3,4-Methylenedioxyamphetamine HCl (MDMA)	50
(±) 3,4-Methylenedioxyamphetamine HCl (MDA)	300
3,4-Methylenedioxyethyl-amphetamine (MDE)	30
l-Methamphetamine	25,000
METHAMPHETAMINE (METS0)	
d-Methamphetamine	50
Fenfluramine	60,000
p-Hydroxymethamphetamine	400
Methoxyphenamine	25,000

Mephentermine	1,500
3,4-Methylenedioxyamphetamine (MDMA)	50
(1R, 2S) - (-)-Ephedrine	400
Procaine	2,000
l-Phenylephrine (R)-(-)-Phenylephrine	6,250
Ephedrine	400
Benzphetamine	25,000
OPIATES (MOR/OPI50)	
Morphine	50
Codeine	65
Ethylmorphine	65
Hydromorphone	250
Hydrocodone	150
Levorphanol	1,000
Oxycodone	75,000
Morphine 3-β-D-Glucuronide	90
Normorphine	90,000
Nalorphine	>100,000
Oxymorphone	65,000
Thebaine	35,000
Diacetylmorphine (Heroin)	65
6-Monoacetylmorphine	150
OPIATES (MOR/OPI40)	
Morphine	40
Codeine	50
Ethylmorphine	50
Hydromorphone	200
Hydrocodone	100
Levorphanol	800
Oxycodone	60,000
Morphine 3-β-D-Glucuronide	70
Normorphine	70,000
Nalorphine	10,000
Oxymorphone	50,000
Thebaine	25,000
Diacetylmorphine (Heroin)	50
6-Monoacetylmorphine	125
OPIATES (MOR/OPI10)	
Morphine	10
Codeine	5
Ethylmorphine	25
Hydromorphone	70
Hydrocodone	270
Levorphanol	1,000
Oxycodone	>10,000
Morphine 3-β-D-Glucuronide	20
Normorphine	10,000
Nalorphine	700
Oxymorphone	>10,000
Thebaine	>10,000
Diacetylmorphine (Heroin)	25
6-Monoacetylmorphine	10
METHADONE (MTD30)	
Methadone	30
Disopyramide	400
(+)-Chlorpheniramine	6,250
LAAM	200
Doxylamine	12,500
Nor-LAAM	12,500
OXYCODONE (OXY40)	
Oxycodone	40
Oximorphone	80
Levorfanol	20,000
Hydrocodone	3,000
Hydromorphone	20,000
Codeine	50,000
Dihydrocodeine	12,500
Naloxona	10,000
Naltrexona	10,000
Tebaina	50,000
OXYCODONE (OXY20)	
Oxycodone	20
Oximorphone	40
Levorfanol	10,000
Hydrocodone	1,500
Hydromorphone	10,000
Codeine	25,000
Dihydrocodeine	6,250
Naloxona	5,000

Naltrexona	5,000
Tebaina	25,000
PHENCYCLIDINE (PCP10)	
Phencyclidine	10
4-Hydroxyphencyclidine	2,500
SYNTHETIC CANNABINOIDS (SPC25)	
JWH-018 5-Pentanoic acid	25
JWH-073 4-Butanoic acid	25
JWH-018 4-Hydroxypentyl	210
JWH-018 5-Hydroxypentyl	300
JWH-073 4-Hydroxybutyl	170
JWH-018 N-Propanoic acid	20
JWH-019 6-Hydroxyhexyl	500
JWH-122 N-4-Hydroxypentyl	500
RCS4 N-5-Carboxypentyl	22,500
MAM2201 N-Pentanoic acid	35
JWH-210 N-5-Carboxypentyl	210
JWH-398 N-Pentanoic acid	175
JWH-200 6-Hydroxyindole	300
JWH-073 N-2-Hydroxybutyl	500
JWH-019 5-Hydroxyhexyl	500
JWH-018	42,000
AM2201 N-(4-hydroxypentyl)	350
JWH-073 N-(3-hydroxybutyl)	225
CANNABIS (THC40)	
Δ ⁹ -THC	40
Cannabinol	40,000
(±)-11-Hydroxy-Δ ⁹ -THC	800
11- nor-Δ ⁹ -THC-9 COOH	32
(-) Δ ⁸ -THC	250
(±) Δ ⁸ -THC	80
CANNABIS (THC15)	
Δ ⁸ -THC	15
Cannabinol	20,000
(±)-11-Hydroxy-Δ ⁹ -THC	400
11- nor-Δ ⁹ -THC-9 COOH	12.5
(-) Δ ⁸ -THC	100
(±) Δ ⁸ -THC	40
TRAMADOL (TML30)	
Cis-tramadol	30
Procyclidine	3,000
d,l-O-Desmethyl venlafaxine	15,000
n-Desmethyl-cis-tramadol	15
Phencyclidine	6,000
o-Desmethyl-cis-tramadol	1,500

(+/-) Brompheniramine	Ethanol (Except ALC)	Norethindrone	Thebaine
Benzocaine	Etodolac	Noscapine	Theophylline
Buspirone	Famprofazone	Niacinamide	Thiamine
Caffeine	Fenoprofen	Norephedrine	Thioridazine
Chloramphenicol	Fluoxetine Hydrochloride	Orphenadrine	Tolbutamide
Chloroquine	Furosemide	Oxalic acid	Trazodone
(+/-)-Chlorpheniramine	Gentisic acid	Oxolinic acid	Triamterene
S- (+)-Chlorpheniramine maleate salt	D (+) Glucose	Oxymetazoline	Trifluoperazine
Chlorpromazine	Guaiacol Glyceryl Ether	Papaverine	Trimethoprim
Chlorprothixene	Hemoglobin	Pemoline	Trimipramine
Cimetidine	Hydralazine	Penicillin-G	Tryptamine
Clomipramine	Hydrochlorothiazide	Perphenazine	Tyramine
Clonidine	Hydroxyzine	Phenelzine	Uric acid
Creatine	Imipramine	Pheniramine	Verapamil
Cyclobenzaprine	Isoproterenol hydrochloride	Phenothiazine	Zomepirac

Literature

1. Moolchan, E., et al, "Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine", Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.

2. Kim, I, et al, "Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration", *ClinChem*, 2002 Sept.; 48 (9), pp 1486-96.

3. Schramm, W, et al, "Drugs of Abuse in Saliva: A Review," *J Anal Tox*, 1992 Jan-Feb; 16 (1), pp 1-9

4. McCarron, MM, et al, "Detection of Phencyclidine Usage by Radioimmunoassay of Saliva," *J Anal Tox*. 1984 Sep-Oct.; 8 (5), pp 197-201.

Rev.1.03; 2025-01-22 MaKa/BIrU

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds spiked into drug-free PBS stock. The following compounds demonstrated no false positive results on the nal von minden Drug-Screen® Saliva Tracer when tested with at concentrations up to 100 µg/mL.

Acetaminophen	Dextromethorphan	Isoxsuprine	β-Phenylethylamine
Acetone	Diclofenac	Kanamycin	Procaine
Acetophenetidin	Dicyclomine	Ketoprofen	Promethazine
Aspirin	Diffunisal	Labetalol	Quinacrine
Albumin	Digoxin	Lidocaine	Quinidine
Amoxapine	4-Dimethylaminantipyrine	Lindane	Ranitidine
Amoxicillin	Diphenhydramine	Loperamide	Riboflavin
Ampicillin	5,5-Diphenylhydantoin	Meperidine	Sodium chloride
Ascorbic acid	Disopyramide	Methoxyphenamine	Sulfamethazine
Aspartame	Doxylamine	Metoprolol	Sulindac
Atropine	Dopamine	Nalidixic acid	Temazepam
Benzoic acid	(1R, 2S) - (-)- Ephedrine	(+)-Naproxen	Tetracycline
Bilirubin	Erythromycin	Nimesulide	Tetrahydrozoline

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń
	Bei beschädigter Verpackung nicht verwenden	Do not use if package is damaged	Ne pas utiliser si l'emballage est endommagé	No utilizar si el envase está dañado	Non utilizzare se la confezione è danneggiata	Nie używać w przypadku uszkodzonego opakowania
	Bevollmächtigter in der Europäischen Gemeinschaft	Authorised representative in the European Community	Mandataire dans l'Union Européenne	Representante autorizado en la Comunidad Europea	Mandatario nella Comunità Europea	Upoważniony przedstawiciel we Wspólnocie Europejskiej
	Schweizer Bevollmächtigter		Mandataire suisse		Mandatario svizzero	
	Distributor	Distributor	Distributeur	Distribuidor	Distributore	Dystrybutor

Symbol	Português	Český	Suomi	Svenskt	Dansk	Nederlands	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-mærkning	CE-markering	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyt-töohjetta	Läs bruksanvisningen	Se brugsanvisnin-gen	Raadpleeg de ge-bruiksaanwijzing	Les bruksanvis-ning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdra-votnický pros-tředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	<i>in-vitro</i> diagnostic medisk enhet
	Limites de tempe-ratura	Teplotní omezení	Lämpötilarajat	Temperaturbe-gränsning	Temperaturbe-grænsning	Temperatuurlimiet	Temperatur be-grensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Batchkode	Code van de partij	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvån-das	Må ikke genbruges	Niet opnieuw ge-bruiken	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Udløbsdato	Houdbaar tot	Tidtaking
	Número de catálogo	Katalogov éčíslo	Lueteltonumero	Listnummer	Best il l ingsnum-mer	Catalogus nummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Tilstrækkeligt til <n> test	Voldoende voor <n> test	Tilstrækkelig for<n> tester
	Não utilizar se a embalagem estiver danificada	Nepoužívejte, je-li obal poškozen	Käyttö kielletty jos pakkaus on vaurioitunut	Får inte användas om förpackningen är skadad	Må ikke anvendes, hvis emballagen er beskadiget	Niet gebruiken indien de verpakking beschadigd is	Skal ikke brukes hvis pakningen er skadet
	Mandatário na Co-munidade Euro-peia	Zplnomocněný zástupce v Evropském společenství	Euroopan yhteisössä toimiva valtuutettu edustaja	Auktoriserad representant inom Europeiska gemenskapen	Angiver den autoriserede repræsentant i Den Europæiske Union	Gemachtigde in de Europese Gemeenschap	Autorisert representant i EU
	Distribuidor	Distributor	Jakelija	Distributör	Distributør	Distributeur	Distributør

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info@nal-vonminden.com

www.nal-vonminden.com

Hangzhou Biotest Biotech Co., Ltd.
17#, Futai Road, Zhongtai Street,
Yuhang District, Hangzhou, P.R. China

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80,
20537 Hamburg, Germany

SUNGO Technical Service GmbH
Bahnhofstrasse 21,
6300 Zug, Switzerland

nal von minden GmbH
Carl-Zeiss-Strasse 12 • 47445 Moers • Germany
www.nal-vonminden.com • info@nal-vonminden.com