nal von minden Drug-Screen[®] Kratom (KRA) Test (Urine)

REF. 101099, 102099

DE	Produktinformation und Gebrauchsanleitung	2	FI	Tuotetieto ja käyttöohje	18
EN	Product information and instructions for use	4	SE	Produktinformation och bruksanvisning	20
FR	Informations sur le produit et instructions d'utilisation	6	DK	Produktinformation og brugervejledningen	22
ES	Información del producto e instrucciones de uso	8	NL	Productinformatie en instructies voor gebruik	24
IT	Informazioni ed istruzioni per l'uso	10	NO	Produktinformasjon og bruksanvisning	26
PL	Informacje o produkcie i Sposób użycia	12		Performance Characteristics	28
PT	Informação sobre o produto e instruções de utilização	14		Literature	28
CZ	Informace o produktu a návod k použití	16		Symbols	31
				Our Teams	32

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1. Intended Use

The nal von minden Drug-Screen® Kratom (KRA) Test is an immunological rapid test intended for the qualitative determination of mitragynine in human urine:

Parameter*	Drug/Metabolite**	Cut-off [ng/mL]
KRA	Mitragynine	100

*Abbreviation printed on test.

**Substance used for calibrating the different cut-off levels.

The nal von minden Drug-Screen® rapid tests are *in-vitro* diagnostic devices and are intended for professional use only. The rapid tests only provide a visual, qualitative and preliminary analytical test result. A secondary analytical method should be performed to confirm results. Gas chromatography-mass spectrometry (GC-MS) or liquid chromatography-mass spectrometry (LC-MS) are the preferred confirmatory methods. Clinical considerations and a professional evaluation of every test result should be obtained, especially in cases of preliminary positive results.

2. Test Principle

All nal von minden Drug-Screen® drug rapid tests are based on the same, approved test principle. The patient sample is applied to the sample pad of the test and moves along the test strip by capillary action. In this way, free, gold-conjugated antibodies, located near the sample area, are transported along the test strip along with the urine sample. In the event of a negative sample, these reach the test line region, where the test's target substance (e. g. THC) is immobilised. These immobilised drugs are recognised and bound by the gold-conjugated antibodies. As a result, gold particles accumulate at the test line region (T) causing the formation of a red line, which indicates a negative test result. If a drug is present in the sample, it is bound after the application of the sample by the gold-conjugated antibodies in the conjugate pad. If the drug concentration in the patient sample is higher than the cut-off level of the test strip, all binding sites of the gold-conjugated antibodies are saturated and cannot be bound at the T-line. In the event of a positive sample above the cut-off level, the T-Line fails to appear. Regardless of the presence of a drug, gold-conjugates will reach the control line region and will be bound there by the antibodies on the membrane. To indicate that a test is valid, a red control line will always develop.

3. Material

Provided:

- Single tests (Test cassettes incl. disposable pipettes)
- Package insert

Additional materials required:

- Timer
- Gloves
- Urine collection cup

4. Storage and Stability

nal von minden Drug-Screen® rapid tests can be stored either at room temperature or refrigerated (2-30°C) at normal humidity. Do not use the tests beyond the stated expiry date. The product is sensitive to humidity. Individually packed tests should be used immediately after opening the individual

pouch. Tests from damaged pouches should be discarded. When single test strips packed in a tube are used, the tube should be resealed securely after the removal of a test strip. Tests must be used within 3 months of opening the tube. Please make a note of the date of opening.

5. Warnings and Precautions

- For *in-vitro* diagnostic use only.
- For professional use only.
- For single use only.
- Store tests between 2-30°C. Do not freeze.
- Do not use tests beyond the expiry date.
- Do not use tests with damaged pouches.
- Once removed from its foil pouch, the test should be used immediately.
- Do not touch the sample area or the reaction fields/pads.
- Sample materials are potentially infectious. During the test procedure, adhere to standard guidelines for the handling of potentially infectious material and chemical reagents. The use of protective clothing (laboratory coat, gloves, eye protection) is recommended. Materials that have come into contact with the sample material should be disposed of in accordance with local regulations.

6. Specimen Collection, Preparation and Storage

The urine specimen must be collected in a clean and dry plastic or glass container. Urine collected at any time of day may be used. If testing is not to be carried out immediately, the urine specimens may be stored at 2-8°C for up to 48 hours. For long-term storage, specimens should be kept at or below -20°C.

Refrigerated urine samples should be brought to room temperature prior to testing. Frozen specimens must be mixed thoroughly after thawing. Repeated freeze/thaw cycles should be avoided.

7. Test Procedure

Refrigerated tests and urine samples should be brought to room temperature (15-30°C).

1. Remove the test from the foil pouch.
2. **A) Dip test:**

Immerse the test strip in the urine sample as far as the corresponding mark (**MAX**) for **10-15 seconds**. The liquid must not exceed the MAX mark.

B) Cassette test:

Transfer **3 drops of urine** (about 120 µL) into each sample well of the test cassette. Ensure that the urine does not come into contact with other areas of the test.

A) Dip test

B) Cassette test

3. Place the test on a clean and flat surface.

4. Start the timer.
5. Reading the results:
 - a. Read the results of the drug tests after **5 minutes**. **Do not read the results after more than 10 minutes**.
 - b. If using a Rapid Slide Scanner (RSS), interpret the tests at the indicated time. For this purpose, you can also use the integrated timer function. Ensure you select the correct test using the scanner-software (see RSS manual).

- Test results should always be interpreted in the light of all the evidence and never in isolation.
- Results determined as positive indicate the presence of the corresponding drugs/medications in urine, but do not verify poisoning or its extent, and they cannot comment on the frequency or level of consumption.
- There is a possibility that the test result is distorted by technical error, error in the test procedure or by substances or factors that influence the test and have not been mentioned here.

8. Result Interpretation

The reaction field contains a test line region (T) for the drugs that are to be detected, as well as a control line region (C).

Rev.1.00 2022-05-23 HeMu

Negative


A drug-free urine sample or a urine sample with drug content below the detection limit (cut-off) of the test was analysed. A test line appears in test line region (T).

Positive

No line appears in the test line region (T).

Invalid

No control line appears in the control line region (C). The test result should be disregarded. The test must be repeated with a new test device.

Note:

The colour intensity of the C- and T-lines on the test membrane may vary. Only if no T-line appears is the test positive. Positive and unclear results should be confirmed with another analytical method (e.g. GC-MS).

9. Quality Control

With the control line (C), the nal von minden Drug-Screen® drug rapid tests have an integrated process control. The control line is formed as a result of an independent antigen/antibody reaction and must always appear independent of drug and metabolite concentrations in the sample. The control line confirms that sufficient sample material has been added and the test was performed correctly. For this reason, we recommend recording the appearance of the control line when performing the test. If no control line appears, the test is invalid and must be disposed of. Please review the procedure and repeat the test with a new test. If the problem persists, discontinue using the batch immediately and contact your distributor.

10. Test Limitations

- The nal von minden Drug-Screen® drug tests are only suitable for the analysis of human urine.
- Cross-reactivity or interaction profiles must be considered in the evaluation of the nal von minden Drug-Screen® drug tests. Please refer to the 'Analytical Specificity' section at the end of the package insert.
- Positive results obtained with the nal von minden Drug-Screen® drug tests must be confirmed using another analytical method.

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Performance Characteristics

Analytical Specificity

Drug – Cut-off 100 ng/mL	Concentration (ng/mL)
Mitragynine	100
7-Hydroxymitragynine	125
Olanzapine	50,000

With exception of the respective parameters positive reacting drugs and drug metabolites, which are listed above, all below listed compounds react negatively up to a concentration of 100 µg/mL.

Acetaminophen	4-Dimethylamino-antipyrine	Niacinamide
Acetone	Diphenhydramine	(+/-)-Norephedrine
Albumin	Dopamine	Oxalic acid
Ampicillin	(+/-)-Isoproterenol	Penicillin-G
Ascorbic acid	(+)-Naproxen	Pheniramine
Aspartame	Erythromycin	Phenothiazine
Aspirin	Ethanol	L-Phenylephrine
Atropine	Furosemide	B-Phenylethylamine
Benzocaine	Glucose	Procaine
Billirubin	Gualacol glyceryl ether	Quinidine
Caffeine	Hemoglobin	Ranitidine
Chloroquine	Ibuprofen	Riboflavin
(+)-Chlorpheniramine	(+/-)-Isoproterenol	Sodium chloride
(+/-)-Chlorpheniramine	Levorphanol	Sulindac
Creatine	Lidocaine	Theophylline
Dexbrompheniramine	(1R,2S)-(-)-n-Methylephedrine	Tyramine
Dextromethorphan		

Analytical Sensitivity

Drug-free urine and urine with drug concentrations at +/-50% cut-off (50 ng/mL and 150 ng/mL) and +/-25% cut-off (75 ng/mL and 125 ng/mL) were tested by 9 operators at 3 physician office laboratories over 20 non-consecutive days or by in-house personnel at the manufacturing site. Results showed over 99% agreement at +/-50% cut-off levels with the nal von minden Drug-Screen® KRA100 Test.

Rev.1.00 2022-05-23 UISi

Literature

- Stephanie Colucci. Why Dose the Food and Drug Administration Need to Ban Kratom[J]. Current Emergency and Hospital Medicine Reports, 2019, Vol.7 (4), pp.169-174.
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488.
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

Rev.1.00 2022-05-23 UISi

Specimen Correlation

The accuracy of the nal von minden Drug-Screen® KRA100 Test was evaluated in comparison to GC/MS. Drug-free urine samples collected from presumed non-user volunteers were tested with KRA Rapid Test. Of these negative samples, all were correctly identified as negative. 10% of the negative samples were confirmed with GC/MS as drug negative. 21 positive urine specimens were obtained from reference labs. Drug concentrations were confirmed with GC/MS. A summary of the accuracy results are summarized below:

Drug-Screen KRA100 Test	result	Range of GC/MS Data					
		Drug free	50% - 75% C/O	75% - 100% C/O	100% - 125% C/O	125% - 150% C/O	> 150% C/O
Test strip	Neg.	40	2	0	0	0	0
	Pos.	0	0	1	1	3	14

% Positive agreement: 97.67%

% Negative agreement: >99%

% Total agreement: 98.36%


Drug-Screen KRA100 Test	result	Range of GC/MS Data					
		Drug free	50% - 75% C/O	75% - 100% C/O	100% - 125% C/O	125% - 150% C/O	> 150% C/O
Test cassette	Neg.	40	2	0	0	0	0
	Pos.	0	0	1	1	3	14

% Positive agreement: 97.67%

% Negative agreement: >99%

% Total agreement: 98.36%

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limitación de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour pour "n" tests	Válido para para <n> ensayos	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återvändas	Niet opnieuw gebruiken	Må ikke genbruges
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato
	Número de catálogo	Katalogov číslo	Luettonumero	Listnummer	Catalogus nummer	Bestilingsnummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax 0800 024 9519

Denmark:

Tel: +31 703075 603
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH