

PEREG GMBH
PORSCHESTRASSE 12
D-84478 WALDKRAIBURG
TELEFON (0 86 38) 8 41 95
TLEFAX (0 86 38) 8 41 62
E-MAIL pereg@t-online.de
http://www.pereg.de

Declaration of conformity EN ISO 15883-5 2021

PEREG GmbH protein residue test:

Pyromol-Test and Pyromol-E

EN ISO 15883-5 outlines examples of test methods for the detection of proteins (ANNEX C).

Other protein test methods can be used provided they are validated.

The validated Pyromol-Test and Pyromol-E are based on the formation of a dye protein complex and can detect insoluble protein residue down to 1µg.

The standardised PEREG residue tests are produced to enable swabbing the targeted areas. The detection level is below the alert level assay criteria outlined in 4.4.3 and the qualitative result can be used for routine testing.

PEREG GmbH

Martin Pfeifer
Managing Director

PEREG GMBH

PORSCHESTRASSE 12 D-84478 WALDKRAIBURG TELEFON (0 86 38) 8 41 00 TELEFAX (0 86 38) 8 41 62

January 2022