

# Rapid Slide Scanner (RSS) III – automatic evaluation and documentation of rapid tests and results

The RSS III is the ideal tool for the objective interpretation of all kinds of rapid tests. It evaluates not only test strips, but also test cassettes and test cups within a matter of seconds, delivering clear results and comprehensive documentation.

The RSS III software recognises the position and format of the test, so that it can be interpreted wherever it is placed on the scanner. Test cassettes with barcodes are automatically identified and evaluated.

## Advantages

- ✓ Tests can be positioned anywhere on the scanner
- ✓ Objective evaluation in a matter of seconds
- ✓ Automatic barcode recognition
- ✓ Software in multiple languages
- ✓ Simultaneous evaluation of up to 12 parameters
- ✓ Individual result reports
- ✓ Reliable, competent customer support
- ✓ LED technology
- ✓ Remote maintenance and updates via internet
- ✓ Barcode scanner (handheld; optional) for the acquisition of patient data
- ✓ LIS connection (optional)
- ✓ Data export via Microsoft-Excel™ or CSV (optional)
- ✓ Data import via CSV (optional)

